

opSens<sup>...</sup>



# AccuSens

## User Manual



# AccuSens User manual

## Opsens inc

2014 Cyrille Duquet, Suite 125  
Sainte-Foy, Quebec  
G1N 4N6, Canada  
Tel. : +418-682-9996  
Fax: +418-682-9939  
[info@opsens.com](mailto:info@opsens.com)  
[www.opsens.com](http://www.opsens.com)

## **WARRANTY**

All products manufactured by Opsens inc. are warranted to be free of defects in workmanship and materials for a period of one year from the date of shipment. No other express warranty is given, and no affirmation of Seller, by words or actions, shall constitute a warranty. SELLER DISCLAIMS ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSES WHATSOEVER. If any defect in workmanship or material should develop under conditions of normal use and service within the warranty period, repairs will be made at no charge to the original purchaser, upon delivery of the product(s) to the factory, shipping charges prepaid. If inspection by Opsens or its authorized representative reveals that the product was damaged by accident, alteration, misuse, abuse, faulty installation or other causes beyond the control of Opsens, this warranty does not apply. Service, repairs or disassembly of the product in any manner, performed without specific factory permission, voids this warranty.

OPSENS MAKES NO WARRANTY OF ANY KIND WITH REGARD TO THIS MANUAL, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. Opsens shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

This warranty does not apply to the transducers sold for use with Opsens' signal conditioners.

## **IMPORTANT NOTICE**

The product specifications and other information contained in this manual are subject to change without notice.

Opsens has made a concerted effort to provide complete and current information for the proper use of the equipment. If there are questions regarding this manual or the proper use of the equipment, contact Opsens inc at:

TEL + 418-682-9996  
FAX + 418-682-9939

WEB SITE www.opsens.com  
E-MAIL info@opsens.com

## **PACKAGING FOR SHIPMENT**

If the product must be shipped to a different location or returned to Opsens for any reason through a common carrier it must be properly packaged to minimize risks of damage.

## **SOFTWARE LICENSE AGREEMENT**

This product contains intellectual property, i.e. software programs, that are licensed for use by the end user/customer (hereinafter "end user").

This is not a sale of such intellectual property.

The end user shall not copy, disassemble or reverse compile the software program.

The software programs are provided to the end user "as is" without warranty of any kind, either express or implied, including, but not limited to, warranties of merchantability and fitness for a particular purpose. The entire risk of the quality and performance of the software program is with the end user.

## WARNINGS



### **WARNING: HIGH PRESSURE!**

High pressure gases and liquids are potentially hazardous. Energy stored in these gases or liquids can be released suddenly and with extreme force. High pressure systems should be assembled and operated only by personnel who have been trained in proper safety practices.



### **WARNING : HIGH STRAIN!**

Highly strained materials and parts are potentially hazardous. Energy stored in these materials or parts can be released suddenly and with extreme force. Highly strained systems should be assembled and operated only by personnel who have been trained in proper safety practices.



### **WARNING: NOT EXPLOSION PROOF!**

Installation of this instrument in an area requiring devices rated as intrinsically safe is not recommended.



### **WARNING: VOLTAGE SUPPLY!**

Use only the wall plug-in power supply delivered with your PicoSens and verifies that the input voltage and frequency are compatible with the power outlet.



## Table of content

<b>Table of content</b> .....	<b>v</b>
<b>1. AccuSens Operating Principle</b> .....	<b>1</b>
<b>2. Quick start</b> .....	<b>2</b>
<b>2.1 Power supply</b> .....	<b>2</b>
<b>2.2 Sensor connection</b> .....	<b>2</b>
<b>2.3 AccuSens setup</b> .....	<b>4</b>
2.3.1 Define a sensor.....	4
2.3.2 Select a sensor among previously defined sensors.....	6
<b>2.4 AccuSens remote connection and control</b> .....	<b>6</b>
<b>2.5 Analog output</b> .....	<b>6</b>
<b>3. Local Operation</b> .....	<b>8</b>
<b>3.1 Keyboard</b> .....	<b>8</b>
3.1.1 Switch on/off .....	8
3.1.2 Menu Button .....	8
3.1.3 Left/Right Arrows .....	8
3.1.4 Up/Down Arrows .....	8
3.1.5 Confirmation button.....	8
3.1.6 Scan Button .....	8
3.1.7 Define Button .....	9
3.1.8 Select Button .....	9
3.1.9 Null Button .....	9
<b>3.2 Main Menu</b> .....	<b>9</b>
<b>3.3 Probes Submenu</b> .....	<b>9</b>
3.3.1 Def Probe.....	9
3.3.2 Sel Probe.....	10
3.3.3 Mod Probe.....	10
3.3.4 Del Probe.....	10
<b>3.4 Setting Submenu</b> .....	<b>10</b>
3.4.1 Channels .....	10
3.4.2 Average .....	10
3.4.3 Analog output.....	10
3.4.4 Units .....	11
3.4.5 Special.....	11
<b>3.5 Diagnostic Submenu</b> .....	<b>11</b>
<b>4. Remote operation</b> .....	<b>13</b>





## 1. AccuSens Operating Principle

The AccuSens is a fiber-optic white light interferometric multi-channel signal conditioner having the capability of accurately measuring the absolute path length difference of various type of sensing interferometers. Opsens produces a variety of interferometric sensors that are based on either so-called Fabry-Perot interferometer (low-finesse version) configuration or on the polarization interferometer configuration. For example, Opsens temperature sensors are based on a polarization sensing interferometer using the temperature-dependent optical properties of a birefringent crystal while Opsens pressure, strain and displacement sensors are based on a Fabry-Perot sensing interferometer where the distance between the two mirrors of the interferometer varies as a function of the measured parameter. In all cases, the sensing interferometer is made so its path length difference varies with the physical parameter of interest. The path length difference of the sensing interferometer is accurately measured with a nanometer resolution and this over 30 000 nanometers range.

Physical measurements are possible if the path length difference of the sensing interferometer is a univocal function of the parameter under scrutiny. This being the case with Opsens sensors, the AccuSens must know the relation between the path length difference and the physical parameter. The Gauge Factors, or equivalently the calibration factors, that comes with Opsens sensor contain all the information needed by the AccuSens to perform the conversion from the measurement of the path length difference to the physical value being measured.

For readers that may need to have a better understanding of the operating principle of the AccuSens, we refer to the white paper IMP0002 entitled “Opsens White Light Polarization Interferometry Technology” available on Opsens’ web site: <http://www.opsens.com/pdf/WLPI.pdf>.

## 2. Quick start

### 2.1 Power supply

The AccuSens accepts supply voltage form 12 to 30 Vdc. The supply connection is made through the round 2.1 mm power jack (outside ground, inside positive) located on the back panel (Figure 1). You can use either the included 12V wall-mount power supply or use another power-supply.




**Figure 1 :** AccuSens backpanel

### 2.2 Sensor connection

Opsens fiber-optic sensors or transducers must be mated to the AccuSens output connectors (Figure 2). The optical connectors provided with the AccuSens is usually a square push-pull SC-type connector mounted with a SC-SC type mating. Remove the protective cap of the mating, engage the sensor connector with the orientation key properly oriented and push until it clicks into its place. Always clean the sensor connector prior to connect it to the signal conditioner with the included connector cleaner. Use a fresh section of cleaning cloth every time. Ask for Opsens Fiber-Optic Cleaning Guide for further information on how to clean fiber-optic connectors.

**NOTE:** Always clean the connector ferrule end face with the include connector cleaner before each connection. Failure to do so can result in permanent damage to the connector caused by hard particles trapped between the fiber optic end faces.

**NOTE:** Always replace the protective dust cap on the mating when there is no sensor connected and always replace the protective dust cap on the sensor fiber-optic connector when not in used.



**Figure 2 :** Remove the protective cap before connecting



**Figure 3 :** Always clean the end of the connector ferrule before each connection




**Figure 4 :** Engage the sensor connector with the orientation key properly oriented and push until it clicks into place



**Figure 5 :** To disconnect, grab the square connector body and just pull

## 2.3 AccuSens setup



The AccuSens is only compatible with Opsens WLPI-type high accuracy temperature sensors. (Note that the GaAs type PowerTemp temperature sensors are compatible with either the PowerSens or the PicoPowerSens signal conditioners only). To properly use a specific sensor, its corresponding Sensor Type and Gauge Factor(s) must be

entered into the AccuSens non-volatile memory (Sensor definition) and then the defined sensor must be selected (Sensor selection) on the specific channel. These numbers contain the parameterization factors of the sensor that is the sensor type and the sensor calibration parameters. These factors are indicated on a label fixed on the optical cable of the sensor, nearby the optical connector as indicated on the left figure. The label shows the identification number of the sensor (used for record purpose only), the Gauge Factors GF0, GF1, ..., and the Sensor Type.

The type of the sensors depends of the physical parameter to be measured. For the Accusens, only temperature can be measured, using the Tq-type sensor. There are also special type definitions that are the N-type and the X-type, which are not related to any specific sensor. These are needed for internal calibration purpose only and should be not used unless instructions have been given by Opsens to do so.





### 2.3.1 Define a sensor

Before selecting a specific sensor, it must be defined that is its corresponding Sensor Type and Gauge Factor must be stored into the non-volatile memory of the instrument. The following example shows how to define a Tq-type temperature sensor.

<b>DEF</b>	Enter Define Sensor Menu	Def Probe >Type Tq
	Press (✓) key to accept the type and (✓) key again to edit the name of the sensor	<pre> Def Probe &gt;Nam Tq00 OK GF0 0000 GF1 0000 GF2 0000                     </pre>
	Use Default Name or select a name for the sensor using up or down arrow keys (00 to 99 numbers only)	<pre> Def Probe &gt;Nam Tq01 OK GF1 0000 GF2 0000 GF2 0000                     </pre>
	Enter Gauge Factor GF1; use arrow keys and validate the entry with the confirm (✓) key	<pre> Def Probe Nam Tq01 OK &gt;GF1 4235 GF2 0000 GF2 0000                     </pre>
	Repeat above procedure for GF1	<pre> Def Probe Nam Tq01 OK GF1 4235 GF2 1743 GF2 0000                     </pre>
	And also for GF2	<pre> Def Probe Nam Tq01 OK GF1 4235 GF2 1743 GF2 5424                     </pre>
	Select "OK" and press the confirm (✓) key	<pre> Def Probe Nam Tq01 &gt;OK GF1 4235 GF2 1743 GF2 5424                     </pre>
	Quit the menu	

### 2.3.2 Select a sensor among previously defined sensors

Once a sensor is defined, it can be selected for a specific channel as indicated below (example for selecting in channel 2).

	Enter Sensor Select Menu	1 >= None 2 = None 3 = None 4 = None
	Select a channel using the direction arrows and then the confirm (✓) key	1 = None 2 =>None 3 = None 4 = None
	Select sensor using the up/down arrows	1 = None 2 >= Tq01 3 = None 4 = None
	Quit the menu	

**NOTE:** Without a sensor connected, the AccuSens shows the message “NoS” or “No Signal” on its display.


## 2.4 AccuSens remote connection and control

The AccuSens comes with a RS-232 serial communication port. The AccuSens is remotely controlled by having a computer linked to this serial port and using the SoftSens software provided with the signal conditioner.

With RS-232 serial communication, simply connect a standard serial cable (all wires straight through) to the AccuSens DB-9 connector and to a remote computer. Then open SoftSens software as indicated in SoftSens user manual.

## 2.5 Analog output

Sensor readings are available for each channel at the analog output terminals, located on the backpanel of the instrument. See section 3.4.3 for configuring the offset and scaling of the voltage output. Figure 6 shows the terminal assignments for the analog output.



**Figure 6 : Analog output terminals assignment**

## 3. Local Operation

### 3.1 Keyboard



#### 3.1.1 Switch on/off



The instrument is switched on/off with this button.

#### 3.1.2 Menu Button



Menu button gives access to the system menus. Once in the menu, this button brings the user one level higher into the menu hierarchy. Once at the root, the system will exit the menu and goes back to measurement display.

#### 3.1.3 Left/Right Arrows



Left/Right arrows allows navigating 1) within the Gauge Factors of a given sensor being defined, selected, modified or deleted, or 2) it allows moving from one digit to the other when a value is being entered.

#### 3.1.4 Up/Down Arrows



Up/Down arrows allows navigating 1) between menu items of a given hierarchical level, or 2) it allows changing a value being entered.

#### 3.1.5 Confirmation button



Confirmation button permits confirming 1) a new value being entered, 2) confirm the selection of a menu item then moving one hierarchical level lower, 3) refreshing displayed value in the case of diagnostic.

#### 3.1.6 Scan Button





Pressing this button stop the scanning measurement over all channels and reads only a single channel at full speed. The user can select the channel by using the up/down/left/right arrows. Pressing the SCAN button once again resumes to the scanning of all enabled channels.

### 3.1.7 Define Button



This button is a shortcut that brings the user directly to the menu item for defining the Gauge Factor of a new sensor.

### 3.1.8 Select Button



This button is a shortcut that brings the user directly to the menu item for selecting a new sensor being used.

### 3.1.9 Null Button



This button allows to cancel most operations taking place, and return directly to the measurement display. This button is also used for zeroing a sensor when the system is not within the menu (measurement display mode).

## 3.2 Main Menu

When the Menu button is pressed, the display shows the following submenus :

- Probes
- Setting
- Diagnostic

The Probes submenu is used to define, select, modify and delete sensor definitions in the sensor list.

The Setting submenu is used to modify the general system settings such as averaging mode and analog output configuration.

The Diagnostic submenu is used to display some quality indicators on the optical signal from each channel. This is useful for troubleshooting purpose.

## 3.3 Probes Submenu

The Probes submenu is where all sensor-related commands are made.

### 3.3.1 Def Probe

This command is used to define the parameters for a new sensor. It has been described in section 2.3.1.

### 3.3.2 *Sel Probe*

This command is used to select a sensor for each channel. It has been described in section 2.3.2.

### 3.3.3 *Mod Probe*

This command is used to modify the parameters for a sensor that has been already defined. It works basically the same way as the Def Probe function described in section 2.3.2.

### 3.3.4 *Del Probe*

This command is used to delete a sensor from the sensor list. Select the sensor name to delete with the up/down arrows and pressing the confirm (✓) key.

## 3.4 Setting Submenu

The functions below are used to setup specific parameters of the instrument. Some of the settings are not necessary related to a specific sensor.

### 3.4.1 *Channels*

This function is used to enable (ON) or disable (OFF) some or all channels. Use the up/down arrow to select the channel, press the (✓) key, use the up/down arrow to select ON or OFF and press the (✓) key again.

### 3.4.2 *Average*




The user can set ON or OFF the averaging mode (ON by default). When on, the measurements either displayed on the front panel display, output on the analog output or on the RS-232 port is the result of the average of 5 sequential measurements.

This parameter is not saved and is lost when the system is switched off.

### 3.4.3 *Analog output*

The analog output parameters comprise the scale factor and the offset. The scale factor corresponds to the physical unit per Volt (unit/V) output, while the offset corresponds to the physical value at which the user wants the analog output to be at zero volt. For example, one may desire to have 10 °C/V, being offset at 5 °C. The analog output voltage is thus given by:

$$\text{Temperature} = [\text{Voltage output}] \times 10 \text{ }^{\circ}\text{C/V} + 5 \text{ }^{\circ}\text{C}.$$

To change the analog scale and offset for a specific channel, first select the proper channel with the up  and down  buttons, then press the confirmation button .

Then navigate between the Scale and Offset values, again with the up, down and confirmation buttons.

Any new scale factor is saved in the non-volatile memory. This value is used whenever the system is switched on and off.

By default, the offset value is 0. A new offset value is not saved and it is lost when the system is switched off, or when the sensor is de-selected.

### 3.4.4 Units

The user can choose what are the units displayed for each given type of sensor. For example, temperature sensors (Tq) can have their measurements displayed in either Celsius or Fahrenheit. This choice affects the values on the display, on the analog output as well as the values sent by the serial port.

Note that when the units are changed, the analog output scale usually needs to be adjusted.

For almost all sensor types, the user can select between two kinds of physical units for outputting the data readings as indicated in the following table.

Sensor type		Physical units	
Tq	Medical Temperature	Celsius (°C)	Fahrenheit (°F)

NOTE: changing the kind of physical units of a specific sensor type will affect all the defined sensors of that type.

### 3.4.5 Special

This submenu contains two special functions : “Unit Mode” and “No Signal”.




The “unit mode” setting overrides the display in physical units and gives a reading of the “optical path difference” of the sensors instead. This is a specialized command that should be used only when properly instructed by Opsens personnel.

There are certain circumstances (i.e. some remote control programs) where it is not desirable to display a No Signal when a sensor is unconnected or defective. Hence, the “No Signal” message can be enabled or disabled. This choice affects the values on the display, on the analog output as well as the values sent by the serial port.

Both these settings are not saved and are lost when the system is switched off.

## 3.5 Diagnostic Submenu

The user can look through a variety of internal parameters for each optical channel. This can be used to diagnostic potential problems with the system or a specific sensor. Typical diagnostic parameters are shown below. Diagnostic values can be refreshed by

pressing the confirmation button . The channel number can be changed with the up  and down  buttons.

Parameters		Description
Lg	2.6V(volt)	Light level
Ga	1.3(no unit)	Amplifier gain
Lm	47% (%)	Lamp driving level
Ct	18 %(%)	Signal contrast
SNR	485(no unit)	Signal quality

The following table shows diagnostic values with good signal, poor signal, or with a broken sensor (“fault”). A fault condition results in a “No Signal” being displayed.

Parameter	Good signal	Poor Signal	Fault
Lg	> 2.2	< 2.2	—
Ga	< 2.0	> 2.0	—
Lm	< 90 %	> 90 %	—
Ct	> 15 %	< 15 %	—
SNR	> 200	< 200	< 100

**NOTE:** Without a sensor connected, the instrument shows the message “NoS” or “No Signal” on its display.


In the unlikely situation that this message appears while a sensor is connected to the unit, take note of the diagnostic parameters and contact Opsens technical support.

During a No Signal condition, the analog output and the serial ports output constant values as follow:

Output	No Signal condition output value
Analog	0 Volt
RS-232	65 536.0

## 4. Remote operation

Opsens AccuSens signal conditioner comes with a RS-232 serial communication interfaces to allow control with a remote computer. The RS-232 interface settings are indicated on the following figure.



The image shows a 'Settings' dialog box with five rows of configuration options, each with a label and a dropdown menu:

Label	Value
Baud rate:	9600
Data bits:	8
Parity:	None
Stop bits:	1
Flow control:	None

**Figure 12: RS-232 interface setting**

The serial interface remote control commands are based on the standard SCPI syntax (**S**tandard **C**ommands for **P**rogrammable **I**nstrumentation). The user can create its own remote control software using the various SCPI commands available for the AccuSens. But for ease of operation, Opsens provides its own control software, called SoftSens, which gives access to all the functionalities of the instrument. See SoftSens user manual for how to remotely control the AccuSens conditioner.

For those who which to develop their own remote control software, ask for Opsens Serial communication user manual to get all the information about serial interfacing with Opsens signal conditioners.