



Virtual Reality Platform

- Virtual Reality and Psychophysiology
- BIOPAC VR platform
- BIOPAC VR demos
- BIOPAC VR application notes
- BIOPAC VR resources: library, interface, hardware




Virtual Reality & Psychophysiology

- Controlled and replicable experimental setups
- Manipulation of the environment (and avatars) that would be impossible or prohibitively expensive in the real world
- Synchronization of the events from the virtual world with the physiological data record allowing accurate and automated data analysis.
- Biofeedback



Synchronization and data transfer over the network and/or the parallel port



* MP150 for Mac strongly recommended because it supports network data transfer and features extensive Specialized Analysis tools.



BIOPAC VR demos

VR01 Acrophobia

VR04 Fear of Flying

VR02 Public Speaking

VR05 Iowa Gambling Task

VR03 Cue Reactivity

These demos are intended to serve as a tutorial on how to construct a virtual reality experiment.

- It can be fully modified and has been designed in a modular format with extensive comments to allow reuse of parts in other experiments.

Code is written in the Python programming language and extensive support on programming with Python is provided in the software package and user forums.


- Python is considered easier to master than C++ and other lower level programming languages

3D models from the demo can be reused within the VR platform (only).



Acrophobia - *fear of heights*

The participant is moving up on a scaffold elevator in an environment designed to emphasize the perception of height. Movement can be controlled by the experimenter and/or the participant's physiological responses.





Acrophobia - details

Purpose: Expose the participant to different heights in a virtual environment and record the resulting physiological responses. The participant is moving up on a scaffold elevator in an environment designed to emphasize the perception of height. The movement can be controlled by the experimenter and/or participant's physiological reactions.

The present demo is based on the work of Wilhelm et al (2005).

There are seven “floors” for the elevator (different height levels). There are two modes of the experiment depending on what controls the vertical movement of the elevator:

- **Participant-controlled:** Progress to the next level is allowed only if there were no SCR responses in the last 20 seconds and the SCL level has not been rising for the last 20 sec.
- **Experimenter-controlled:** Progress to the next level is initiated by a keystroke from the experimenter.

In the participant-controlled mode SCL data is sent from the data acquisition machine (ACQ) to the virtual reality rendering machine (VR). The responses are analyzed in real-time to determine when the participant should move to the next height level.

In both experiment modes digital marker data is sent from VR to ACQ. Key events in the virtual world (i.e. going to a new height level) are marked in the physiological record to facilitate automatic data analysis.

Data analysis: Since levels of the independent variable (height of the participants) are marked in the physiological record, an automated analysis of the collected physiological data can be performed, identifying SCR responses for each height level.



Acrophobia - details

Bibliography:

Mechanisms of Virtual Reality Exposure Therapy: The Role of the Behavioral Activation and Behavioral Inhibition Systems. Applied Psychophysiology and Biofeedback, Vol. 30, No. 3, September 2005 (C 2005) DOI: 10.1007/s10484-005-6383-1. Frank H. Wilhelm,^{1,5} Monique C. Pfaltz,¹ James J. Gross,² Iris B. Mauss,² Sun I. Kim,³ and Brenda K. Wiederhold⁴

The development of virtual reality therapy (VRT) system for the treatment of acrophobia and therapeutic case. Jang DP, Ku JH, Choi YH, Wiederhold BK, Nam SW, Kim IY, Kim SI. IEEE Trans Inf Technol Biomed. 2002 Sep;6(3):213-7.



Virtual Reality Treatment in Acrophobia: A Comparison with Exposure in Vivo
P.M.G.Emmelkamp, M.Bruynzeel, L.Drost, C.A.P.G.van der Mast Cyberpsychology and Behavior, Vol.4, No.3, June 2001, pp.335-341

Treatment of acrophobia in Virtual Reality; the role of immersion and presence. Merel Krijn, Paul M. G. Emmelkamp, Roeline Biemond, Claudius de Wilde de Ligny, Martijn J. Schuemie and Charles A.P.G. van der Mast Behaviour Research and Therapy, 2004 Feb; 42(2):229-239.



Public speaking

In this social anxiety scenario, the participant delivers a speech in front of an audience. The audience attitude can be varied by the experimenter between different states (i.e. friendly, hostile, and indifferent). The number of avatars and room size can be controlled.





Public speaking - details

Purpose: Set up a virtual world where the participant can give a speech in front of a virtual audience and record the physiological responses of the participant as parameters like audience behavior, room size and audience size are varied.

Description: The participant delivers a speech from behind a podium. The speech text is presented on a monitor on the podium and can be scrolled via a joystick. The room size and number of people (and whether the people are present in the room) are conditions set in advance. During the experiment, the researcher can switch the audience state – i.e. make the audience display boredom, disapproval, etc. The audience behavior is defined using avatar animations and is not limited to the animations included with the demo.

Data analysis: The current state of the audience is marked in the physiological data Record, thus enabling automatic scoring of the data.

Bibliography:

Brief Virtual Reality Therapy for Public Speaking Anxiety Dec 2002, Vol. 5, No. 6 : 543 -550 Speaking anxiety using virtual reality for exposure. *Depression and Anxiety* 2005;22(3):156-8

An Experiment on Public Speaking Anxiety in Response to Three Different Types of Virtual Audience DP Pertaub, M Slater, C Barker - *Presence: Teleoperators & Virtual Environments*, 2002

An experimental study on fear of public speaking using a virtual environment. *Cyberpsychol Behav.* 2006 Oct;9(5):627-33. Slater M, Pertaub DP, Barker C, Clark DM.



Cue reactivity

The participant is exposed to a sequence of rooms along a corridor that contain different stimulation environments. Four neutral environments and four smoking environments are already included.



Cue reactivity - details

Purpose: A controlled stimulus presentation where the participant is exposed to environments with different stimuli and the physiological responses are recorded. For example, in the case of addictions, this paradigm can be used to investigate the relationship between craving, physiological response and the way the stimulus is presented in the environment.

Description: The participant is exposed to a sequence of rooms along a corridor that contain different stimulation environments. Four neutral environments and four smoking environments are already included. The number of rooms as well as the objects in the rooms can be modified. The participant does not actively navigate the environment.

Data analysis: Digital markers in the physiological record indicate when the participant is exposed to what condition, thus allowing for an automated data analysis to be performed.

Bibliography:

Virtual Reality Cue Reactivity Assessment: A Case Study in a Teen Smoker

Authors: Bordnick, Patrick; Traylor, Amy; Graap, Ken; Copp, Hilary; Brooks, Jeremy

Source: Applied Psychophysiology and Biofeedback, Volume 30, Number 3, September 2005, pp. 187-193(7)

Virtual Reality Cue Reactivity Assessment in Cigarette Smokers Patrick S. Bordnick, Ken M. Graap, Hilary Copp, Jeremy Brooks, Mirtha Ferrer, Cyberpsychology & Behavior. Volume 8, Number 5, 2005



Fear of flying

The participant is seated in an airplane and experiences normal flight, turbulence, and landing. Tactile feedback is employed to increase the experience of presence.





Fear of flying - details

Purpose: Expose participants to an airplane environment and study their physiological responses to different aspects of the experience of being in an airplane.

Description: Participants are immersed in a virtual environment where they are seated in an airplane and experience normal flight, turbulence, and a landing. Tactile feedback is employed (a low-frequency driver is placed underneath the chair) to increase the experience of presence. The experimenter can trigger certain events (i.e. landing sequence, turbulence, cabin announcements, etc.)

Data analysis: All events are marked in the physiological record thus facilitating automated data analysis.

Bibliography:

Cognitive behavior therapy for fear of flying: Sustainability of treatment gains after September 11. Anderson, P., Jacobs, C. H., Lindner, G. K., Edwards, S., Zimand, E., Hodges, L. F., & Rothbaum, B. O. *Behavior Therapy* 37 (2006) 91-97

Three-Year Follow-Up for Virtual Reality Exposure for Fear of Flying. Wiederhold, B.K. (2003) Wiederhold, M.D. *CyberPsychology & Behavior: The Impact of the Internet, Multimedia and Virtual Reality on Behavior and Society.* Vol 6(4). pp 441-445.

Virtual Reality Treatment of Flying Phobia. Rosa M. Baños, Cristina Botella, Concepción Perpiñá, Mariano Alcañiz, Jose Antonio Lozano, Jorge Osma, and Myriam Gallardo. *IEEE TRANSACTIONS ON INFORMATION TECHNOLOGY IN BIOMEDICINE*, VOL. 6, NO. 3, SEPTEMBER 2002



Iowa gambling task

In this classic experiment, the participant has to choose between decks of cards with different payoffs. The skin conductance response before and after making a choice can be easily analyzed due to the marking of events from the experiment in the physiological record. In addition, the appearance of the decks, assigned probabilities of winning and losing, can be modified.



Iowa gambling task - details

Purpose: The Iowa Gambling Task template is an example of a hypothesis testing tool. It is modeled after the work of Bechara et al (1994). See below for a description of the task:


Quote from Bechara et al (2005):

“The participants are given four decks of cards, a loan of \$2000 facsimile US bills, and asked to play so as to win the most money. Turning each card carries an immediate reward (\$100 in decks A and B and \$50 in decks C and D). Unpredictably, however, the turning of some cards also carries a penalty (which is large in decks A and B and small in decks C and D). Playing mostly from decks A and B leads to an overall loss. Playing from decks C and D leads to an overall gain. The players cannot predict when a penalty will occur, nor calculate with precision the net gain or loss from each deck. They also do not know how many cards must be turned before the end of the game (the game in fact ends after 100 card selections).“

Taken from: *The Iowa Gambling Task and the somatic marker hypothesis: some questions and answers* A. Bechara, H. Damasio, D. Tranel and A.R. Damasio
TRENDS in Cognitive Sciences Vol.9 No.4 April 2005
Both text and screenshot taken from Bechara et al (2005).



Iowa gambling task - details



Anticipatory SCR levels change as a function of the number of trials experienced and result in an increasing disparity between levels observed prior to selecting good vs. bad decks. Bechara et al have proposed that this change in somatic response occurs even before the participants have adequate conscious knowledge of the situation.



Iowa gambling task - details

With minor customization this demo application can be directly applied for research or teaching purposes. It allows the user to test the following:

1. Are somatic responses different before and after good vs. bad decks?
2. Are somatic responses different for more vs. less predictable decks (defined as low vs. high variance in outcome)?
3. How do 1, 2 change as a function of time?

Bibliography:

The Iowa Gambling Task and the somatic marker hypothesis: some questions and answers A. Bechara, H. Damasio, D. Tranel and A.R. Damasio. TRENDS in Cognitive Sciences Vol.9 No.4 April 2005

Insensitivity to future consequences following damage to human prefrontal cortex. Bechara, A., Damasio, A. R., Damasio, H. & Anderson, S. W. (1994) Cognition 50, 7–15.

Do somatic markers mediate decisions on the gambling task? Tomb, I., Hauser, M., Deldin, P. & Caramazza, A. (2002) Nat. Neuroscience. 5, 1103–1104. Cleeremans, A. (2001) in International Encyclopedia of the Social & Behavioral Sciences, eds. Smelser, N. J. & Baltes, P. B. (Elsevier, London), Vol. 4, pp. 2584–2589.

Implicit Learning: News From the Front. Cleeremans, A., Destrebecqz, A. & Boyer, M. (1998) Trends Cogn. Sci. 2, 406–416.



VR Resources (in progress)

- 3D model library for the BIOPAC VR platform.
Models can be used within the VR platform only.

Alcohol	Cardboard box	Pit room
– beer bottle	Chair	Podium
– hard liquor bottle	Classroom	Road
– wine bottle	Computer	– curved segment
Auditorium	Desk	– straight segment
Bathroom	Dice	Sky dome
Bookshelf	Dresser	Soccer ball
Buildings	Living room	Soda can
Bus	Monitor	Sofa
Car	Office Cubicle	Table
Card deck (52 objects)	Overpass Scene	TV

- Sound library



VR Interface



Use the WorldViz Vizard Virtual Reality Toolkit (VRT100W) with your BIOPAC data acquisition and analysis system to synchronize events from the virtual world with the physiological data record, allowing accurate and automated data analysis.

Vizard VR Toolkit is everything you need to build complete, interactive 3D content. Designed for rapid prototyping, Vizard gets you creating fast and provides the resources to deploy even the most challenging applications. With Vizard, even someone with no programming experience can leap into the world of interactive 3D content and soon discover what it's like to have an untethered imagination.

One-year support/
maintenance packages
(including upgrades) available:

- **VRT100W-SUP** – 1 seat
- **VRT100W5-SUP** – 5 seats



VR Hardware

- HDS100 Haptic Delivery System




Haptic delivery system for tactile feedback during VR experiments. The system includes an amplifier that connects to a computer sound card and a pair of actuators that vibrate based on the sound from the sound card. Actuators are placed under chair legs or on a platform and deliver vibrations based on the VR environment—i.e., movement of elevators.

- **RXHDS**—Replacement actuators and isolators for HDS100



VR Hardware

- SDS100 Scent Delivery System



Computer controlled (USB) 8-cartridge scent machine that uses compressed air* to project different scents on cue for a predetermined time...followed by a burst of unscented air to clear for the next scent. Scents are triggered from the virtual reality environment. Dispersed scent covers approximately 3-6 meters in front of unit, depending on how many fans are used.

* Requires companion air compressor, SDSAIR or equivalent



VR Hardware

- HMD1 Head-mounted Display



Two high-contrast microdisplays (SVGA 3D OLED) deliver fluid full-motion video in more than 16.7 million colors. The highly responsive head-tracking system provides a full 360-degree angle of view and specially developed optics deliver a bright, crisp image with a nearly 40 deg field of view.



VR Hardware

- HMD2 Head-mounted Display (high res)



State-of-the-art head-mounted display (HMD) for advanced virtual reality applications. Incorporates high-resolution color microdisplays with custom engineered optics to deliver unsurpassed visual acuity in a wide field-of-view format.



VR Hardware

- Precision Position Trackers



Motion capture systems offer high-quality optical tracking over a wide range (more than 10 x 10 m). These PPT systems connect directly to the VR Toolkit for an ideal solution for virtual reality applications.

Using proven CCD technology and

WorldViz image processing technology, the PPT family of products delivers flexible and accurate tracking solutions. Real-time technology displays tracking results the instant a subject performs a motion.

Systems

- **VRPPT2**—2 sensor
- **VRPPT4**—4 sensor
- **VRPPT8**—8 sensor

12-Month Support Packages

- **VRPPT2-SUP**
- **VRPPT4-SUP**
- **VRPPT8-SUP**

Academic discounts available!