



BIOPAC Systems, Inc.

42 Aero Camino
Goleta, CA 93117
Ph (805) 685-0066
Fax (805) 685-0067
www.biopac.com
info@biopac.com

#AS155 - AcqKnowledge File Formats for the Macintosh

Following is the AcqKnowledge file format for the Macintosh platform. This document describes file formatting for AcqKnowledge versions 3.0-3.7.3. This File Format document does not apply to version 3.8.x.

- This document has been formatted for printing and is best viewed as a print out.

```
// BIOPAC Systems AcqKnowledge file format          V3.0                      Nov 20, 1993

// All short doubles are 64 bit IEEE floating point form

// int, & Word are 16 bit, Motorola format (msb first)

// Bool, char & Byte are 8 bits

// Main header

#define nMtMax 6

unsigned int    mainHeaderLen;          // Not used

long           version;                 // File version number:
// 30 = pre-version 2.0, 32 = v 2.0, 33 = v 2.0.7,
35 = v3.0

long           extItemHeaderLen;       // Item header length

int            nChans;                  // Number of channels stored (max = 60)

int            horizAxisType;          // 0 = time, 1 = Freq, 2 = HH:MM:SS, 3=Arbitrary

int            curChannel;              // currently selected channel

short double   sampleTime;              // Sample interval                (msec/sample)

short double   tOffset;                 // time offset                        (msecs)
```

```
short double    tScale;                // time scale                (msecs/div)
short double    tCursor1;              // Cursor 1 time position    (msec)
short double    tCursor2;              // Cursor 2 time position    (msec)

Rect            windRect;               // Window rectangle

int             mmt[nMtMax];            // Measurement functions(1st menu = 1), (pre-
version 3.0 only)

Bool           hilite;                  // Gray non-selected waves
char           dummy1;                  // Not used
short double    firstTOffset;           // Global time offset (msec)
int            rescale;                 // 0=none,1=Autoscale,2=tile; after modify
char           szHorizUnits1[40];       // Horizontal units text (long)
char           szHorizUnits2[10];       // Horizontal units (short)
int            inMemory;                // If non-zero, Keep data file in memory

Bool           grid;                    // Enable grid display
Bool           markers;                  // Enable marker display
int            plotDraft;                // Enable draft mode plotting
Bool           dispMode;                 // 0 = scope mode, 1 = chart mode
char           dummy2;                  // Not used
int            overWritePrompt;         // Not used

// v3.0
int            bShowToolBar;             // Tool bar visible if TRUE
int            bShowChannelButtons;      // Channel numbers visible if TRUE
int            bShowMeasurements;        // Measurements visible if TRUE
int            bShowMarkers;             // Markers visible if TRUE
int            curXChannel;               // Current horizontal axis channel
```

```
int          mmtPrecision;          // Number of digits after decimal point
int          nMeasurementRows;      // Number of measurement rows
int          mmt[40];               // Measurement functions
int          mmtChan[40];           // Measurement channels

// The following structure repeats for nChans as specified above
{
    long      chanHeaderLen;         // Length of channel header
    int       chanNum;               // This channel's number
    char      szComTxt[40];          // Waveform comment text (label)
    long      waveColor;             // waveform color
    int       dispChan;              // 1=invisible, 2 = visible

    short double vOffset;            // ampl offset (units)
    short double vScale;             // ampl scale (units/div)

    char      szUnitsTxt[20];        // ampl units text
    long      bufLength;             // samples
    short double amplScale;          // (units/count) for integer conversion
    short double amplOffset;        // (units) for integer conversion

    int       chanOrder;             // Channel display order (1 is top)
    int       dispSize;              // Channel area display width
    RGBColor  newWaveColor;          // RGB waveform color value
    short     plotMode;              // 0 = scope, 1 = chart, 2 = X-Y
    double    vMid;                  // Midpoint of waveform extents (units)
}
```

```

// End of main header

// Creator specific header

int      creatorHeaderLen;          // Creator header length

int      creatorHeaderType;        // 0x0100 = AcqKnowledge

char     creatorData[creatorHeaderLen-4]; // Creator specific data

// End of creator specific header

typedef struct {
    int     dSize; // Data size (bytes)

    int     dType; // Data type
                // 1 : Floating point
                //           dSize is precision (Bytes)
                //           ( 4, 8, 12, ... )
                // 2 : Binary
                //           dSize is word length

    } binSizeType;

// Data header

binSizeType  dataHead[nChans];      // Data size/types

// Actual data

char         data[]                 // Data is interleaved (lowest number channel first)

{

    long      length;                // Total marker storage length

    long      nMarkers;              // Number of markers stored

```

```
// The following structure repeats for the number of markers as specified above (nMarkers)

{

    long          sample;

    Bool          selected;

    Bool          textLocked;

    Bool          posLocked;

    char          dummy3;          // Not used

    int           textLen;

    char          markerText[textLen];  // Marker text string

}

}
```