

Application Note AH-187 Electrodermal Response (EDR) with an MP System

This application note is concerned with how to record galvanic skin response (GSR), also known as the electrodermal response (EDR), and observe changes with the amplifier module (GSR100C) or the telemetry module (TEL100C) of the MP System (with MP100 or MP150). With the MP System, you will measure both subtle and sudden changes in electrodermal activity, based on your protocols.

Electrodermal Response Basics

Electrodermal response (EDR) measurements show the activity of the eccrine sweat glands. Typically, one will place electrodes where the concentration of these glands is the highest: namely, the fingertips. Responses are a function of the pre-secretory activity of eccrine sweat glands and the filling of the sudorific tubules. The combination of these sudorific elements serves to increase the conductivity of the skin when activated.

When one applies a very small electric voltage (0.5 V) between two electrodes (Ag-AgCl), the manifested electrical conductance varies in direct proportion to the electric current flowing between the electrodes. The electrical conductance is a function of increasing eccrine activity. For instance, if a subject is presented a stimulus and the palms start to sweat, this response indicates a highly-stimulated state. The EDR of this subject will then be higher than the subject's baseline. If another subject receives the same stimulus and the palms remain as "cool as a cucumber," the GSR reading remain unchanged with respect to the baseline. EDR undergoes relatively fast habituation (decrease of amplitude) in the event the same stimulus is repeated over and over to the same subject.

AcqKnowledge software features: EDR

- Automatically locate and analyze data
- Textual event markers
- On-line and off-line analysis
- Remote monitoring through the TEL100C
- On-line journal for note taking
- Automatic stimulation presentation (up to 16 on/off control channels)
- Export results to statistical programs for further analysis
- Determine maximum, minimum, and mean micromho (μmho) measurements

Applications: EDR

- Polygraph ("lie detector")
- Measure stress, arousal, emotional excitement
- Physiological biofeedback
- Measure absolute or relative response levels to different stimuli
- Relaxation Training

MP System Equipment: EDR

EDR data can be collected via the GSR100C or the TEL100C (the wireless model TEL100C-RF is not recommended for EDR applications).

- If using the GSR100C
 - Acquisition Unit: MP100 or MP150
 - Amplifier: GSR100C module
 - Transducer: TSD203 electrodermal response transducer set
 - Electrolyte: GEL101 skin conductance electrode paste or other electrolytic mixture
 - Optional:* Module extension MEC100C (3 meter cable extends lead between GSR100C and TSD203)
- If using the TEL100C
 - Acquisition Unit: MP100 or MP150
 - Amplifier: TEL100C module
 - Transducer: SS3A electrodermal response transducer
 - Electrolyte: GEL101 skin conductance electrode paste or other electrolytic mixture

Hardware Setup

1. Set the Gain switch on the GSR100C.

To measure galvanic skin responses, estimate the approximate units of the skin conductance (μmho) to be investigated.

- The unit of measurement for the GSR100C is the μmho (micromho). Note: $\mu\text{mho} = \mu\text{siemens}$. The mho is the reciprocal of the unit of measurement for resistance, the ohm.
- Larger values indicate a higher level of conductivity; smaller values indicate less conductivity. For example.
 - If the subject is shown provocative slides, the response may be in the 0-50 μmho range.
 - If the subject is presented with a sudden 120dB buzzer in a quiet room, the response may be in the 0-200 μmho range.

Once the appropriate range for the experiment has been determined, choose the proper switch setting from the chart below.

| Conductance Range | | GSR100C Gain Switch |
|-----------------------|-------------------------|----------------------|
| DC | 0.05 Hz | |
| 0-200 μmho | $\pm 200 \mu\text{mho}$ | 20 $\mu\text{mho/V}$ |
| 0-100 μmho | $\pm 100 \mu\text{mho}$ | 10 $\mu\text{mho/V}$ |
| 0-50 μmho | $\pm 50 \mu\text{mho}$ | 5 $\mu\text{mho/V}$ |
| 0-20 μmho | $\pm 20 \mu\text{mho}$ | 2 $\mu\text{mho/V}$ |

For example, for the 0-100 μmho range, set the module Gain switch to the 10 $\mu\text{mho/V}$ setting (on the module Gain switch, mho is an upside down ohm symbol).

If using the TEL100 Remote Monitoring System, for GSR measurements, the following gain settings correspond to μhos . These settings can then be used to calibrate the signal using the rescaling feature in *AcqKnowledge* (see AH-103 for more details on TEL100 setup).

| TEL100 Gain | $\mu\text{mho/V}$ |
|-------------|-------------------|
| 50 | 200 |
| 100 | 100 |
| 200 | 50 |
| 500 | 20 |
| 1,000 | 10 |

| TEL100 Gain | $\mu\text{mho/V}$ |
|-------------|-------------------|
| 2,000 | 5 |
| 5,000 | 2 |
| 10,000 | 1 |
| 20,000 | 0.5 |
| 50,000 | 0.2 |

2. Set the three Filter switches on the GSR100C to the appropriate filtering option for your protocol.

| | | |
|-----------------------|------------|--|
| Low Pass(LP) | 1 Hz LP | Sufficient for almost all EDR studies |
| | 10 Hz LP | Use to investigate higher frequency components of the EDR. |
| DC | DC | Gives direct (absolute) EDR readings from the subject. |
| High Pass (HP) | 0.5 Hz HP | Provides relative EDR recordings. Removes low frequency signals. The effect of using the 0.5 Hz HP setting will be that the subject's baseline will return to nearly zero 1 second after a response change. The result provides an indicator of EDR changes as opposed to a specific EDR level. |
| | 0.05 Hz HP | Provides relative EDR recordings. Removes very low frequency signals. The effect of using the 0.05 Hz HP setting will be that your subject's baseline will return to nearly zero 10 seconds after a response change. The result provides an indicator of EDR changes as opposed to a specific EDR level. |

3. Set the **Channel** switch on the top of the GSR100C.
 - Select 1-16 to correspond with an available Analog Input channel.
4. Connect the **TS203** transducer to the GSR100C.
 - a) Plug the two colored electrode cables into the VIN+ and VIN- inputs.
 - Either blue lead can be connected to either VIN input.
 - b) Plug the black cable into the GND.

GROUNDING When using the GSR100C amplifier with other biopotential amplifiers attached to the same subject, do not attach the ground lead from the biopotential amplifier(s) to the subject. The subject is already appropriately referenced (grounded) to the system via the VIN- attachment to the GSR100C. If a biopotential ground is attached to the subject, then currents sourced from the GSR100C will be split to the biopotential amplifier ground lead, potentially resulting in measurement errors. If biopotential amplifiers do require a ground lead, when used with a GSR100C, then the AC lead (CBL205) can be used in series with the biopotential ground lead.

Optional If using the module extension cable for up to three meters of extra distance between the subject and the MP acquisition unit, plug the **MEC100C** into the GSR100C inputs and then plug the TSD203 into the MEC100C inputs.

Subject Setup using TSD203 Electrode

The TSD203 transducer provides a small constant current between the electrode sites when used with the GSR100C, and the measured resistance between the two electrodes constitutes the electrodermal response. When using the TSD203 to measure electrodermal response, you will need to determine the choice of electrolyte. For effective monitoring of local eccrine activity, use a higher impedance electrolyte with hyposaturated electrolyte concentrations of Cl- (on the order of physiological levels). Many of our researchers have found electrolyte mixtures of 0.05M NaCl are optimal; BIOPAC GEL101 skin conductance electrode paste is a convenient option.



GSR100C shown with MEC100C and TSD203


1. Apply GEL101 to the skin at the point of electrode contact and rub it in. The gel must have a chance to be absorbed and make good contact before recording begins.
2. Make sure the cavities of the TSD203 transducer are clean and fill them with GEL101 (or your own electrolyte mixture).
3. Attach the TSD203 electrode to the subject's fingertips as shown in the picture above, and then wrap the electrode around the fingers with the Velcro straps, applying mild tension.
4. Wait 5 minutes (minimum) before starting to record data.

Software Setup

1. Open the Change Scaling Parameters dialog for the selected GSR100C channel (MP menu > Setup Channels>Scaling).
2. Enter **CAL1 Input Value** as 0.
3. Enter **CAL1 Scale Value** as 0.
4. Enter **CAL2 Input Value** as 1.
5. Enter **CAL2 Scale Value** to match the Gain setting on the GSR100C per the chart below.

| GSR100C Gain Switch | Cal 2 Scale value |
|---------------------|-------------------|
| 20 µmho/V | 20 |
| 10 µmho/V | 10 |
| 5 µmho/V | 5 |
| 2 µmho/V | 2 |

6. Enter the Units as μmho .
7. Click OK to establish the settings and close out of the Scaling Parameters dialog.



Scale settings for GSR100C @ 10 $\mu\text{mho/V}$ Gain

For example:


To use the 0-100 μmho range, set the switch on the GSR100C front panel to 10 $\mu\text{mho/V}$, set the channel on top of the module to the corresponding software channel, go to Scaling, and type in the following values for that channel:

- Cal1: Input: 0, Scale 0
- Cal2: Input 1, Scale 10
- Units: micromhos

Optional Calibration

To verify the Gain setting of the GSR100C:

1. Calibrate *AcqKnowledge* for lower frequency response at DC:



In the Scaling window, set the Input voltages so they map to the "DC" conductance ranges indicated by the sensitivity setting.

- Cal 1: Input 0, Scale 0
- Cal2: Input 1, Scale 5
- Units: micromhos

2. Place both lower frequency response (HP) filters on the GSR100C to **DC**.
3. Set the **Gain** switch on the GSR100C to 5 $\mu\text{mho/V}$.
4. Perform measurement with electrodes disconnected.
5. *AcqKnowledge* should produce a reading of 0 μmho . If slightly off, adjust the ZERO trim pot on the top front of the GSR100C to perfectly zero the reading.
6. Insulate a 100kohm resistor and place it from electrode pad to electrode pad resistor must be insulated from fingers). Perform measurement with electrode-resistor setup.
7. *AcqKnowledge* should produce a reading of 10 μmho .