

Application Note 291: Pulse Amplitude Feedback Control Using BIOPAC Basic Scripting

This application note presents a comprehensive description of how to use the BIOPAC Basic Scripting (BBS), *AcqKnowledge*, and MP hardware to implement simple control loops. Sections include Background, Setup, Algorithm, Results, and Appendices A (channel setting for graph template) and B (code).

1.0 Background

This exercise is intended to show that one can implement feedback control to external devices using MP160/AMI100D and BIOPAC Basic Scripting (BBS). This technique is useful for applications that need/require stability in temperature, flow rate, and pressure.

2.0 Setup

Figure 1 shows the schematic for feedback control. The Flow Measurement signal is processed through CH 8 of the AMI100D then filtered via CH 40 (LPF). The signal from CH 40 is sampled in BBS during OnIdleAcquisition and adjusted as necessary to maintain the desired flowrate. The adjusted signal is routed through OUT0 then to the Flow Source control port (Figure 2a). A bleed valves are used to impact flow and allow the control loop to adjust for it (Figure 2b).




Figure 1 Schematic of Feedback Control of Flowrate




Figure 2a Flow meter signal fed to CH 8; OUT0 fed to flow source control port.




Figure 2b Bleed valves to affect flow.

To correctly applied voltage feedback to the flow source, a mapping of flow rate to voltage is needed. Figure 2 shows that mapping. Also, a regression fit to the mapping points is derived for use in the control macro code.




Figure 3 Mapping of Flow Rate-lpm (abscissa) vs Control Voltage-V (ordinate)

For use in the control code, the regression equation is modified as follows, where F = Flow, V = Voltage:

```

;relation between flowrate and voltage set to stimulator lpm --> voltage
;found through regression equation using discrete data points
V = (-0.0000002*F^3) + (0.0002*F^2) + (0.000006*F) - 0.0115
    
```

3.0 Algorithm

The algorithm is as follows:

The user inputs pulse width in seconds and pulse frequency in hertz which is limited by the duty cycle of 75% referenced to the pulse width entered.

The variable P# is the pulse repetition interval where it is used to sample the waveform for pulse amplitude values.

The variable H#, is that part of the time-domain waveform that is of zero amplitude.

Variable G# is the pulse width in seconds with amplitude V.

OnIdleAcquisition

- during an idle portion of the acquisition, selected the entire graph then extract the max value of the selected section of data

if S# = 0

if (Ticks / 60) - T# > H# ;upon the end of the zero value part of the waveform, set V = a value

Select Wave 8

Edit SelectAll

Set HCursor HCursor2 - P#, HCursor2

X = Max ;sample the waveform for amplitude value

R = X / F ;ratio of measured/desired flow rate

GraphJournal Append STR\$(R, 2)

GraphJournal Append "\r"

if R < 0.1

R = 0.1

V = V / R ;adjust voltage value

MP100 Set Out0, V ;apply voltage value to flow source

S# = 1

T# = Ticks / 60

elseif R >= 0.1

V = V / R ;adjust voltage value

MP100 Set Out0, V ;apply voltage value to flow source

S# = 1

T# = Ticks / 60

endif

endif

elseif S# = 1

if (Ticks / 60) - T# > G# ;upon the end of the non-zero value of the waveform, set V = 0

MP100 Set Out0, 0

S# = 0

T# = Ticks / 60

endif

endif

4.0 Results

Figures 4a, 4b show an acquisition where the **peak pulse** flow rate, specified at 35 LPM, is maintained when leaks are introduced into the system.

After initial settling, a bleed valves are opened to allow leakage of the gas and to cause a drop in the peak pulse flow rate to 24 lpm. The control algorithm rightly corrects for this by adjusting the flow source voltage to produce MORE flow to compensate for the loss.

When the bleed valve is closed, the flow meter senses a transitory increase in peak pulse flow rate to 38 lpm to which the algorithm compensates by adjusting the flow source voltage to produce LESS flow.

Figure 4b shows a truncated view of the active channels as the full setup has 8 analog and 12 calculation channels each at 1000 S/s. The feedback technique, with the proper logic for toggling the flow source is unimpeded by these many channels and the high sample rate.




Figure 4a Pulse amplitude recovery (reference red line on CH 8) on startup, bleed valve open, bleed valve closed. Graph file with 2 channels @ 1000 S/s each. Pulse width = 1.5 sec, pulse freq = 0.3333 Hz




Figure 4b Pulse amplitude recovery (reference red line on CH 8) on startup, bleed valve open, bleed valve closed. Graph file with 20 channels @ 1000 S/s each. Pulse width = 1.5 sec, pulse freq = 0.3333 Hz

Appendix A | Channel Settings for Graph Template

Acquisition Settings

MP160

Serial #: MP160 001F54

Graph: Feedback Flow Ref pulse control test result.acq

Sampling rate: 1000 samples/sec

Record and replace previous data; save to memory.

Acquisition length: 2 hours (7200001 samples)

Smart Amplifier Configuration:

Acquired channels:

Analog Channel 8

Label: Flow Reference

Scaling: 0.00030518 volts -> 0 lpm, 10.0003 volts -> 300 lpm

Units: lpm

Calculation Channel 0

Label: Flow 10Hz LPF - LT

Units: lpm

Filter - IIR

Source channel: Analog Channel 8

Low pass at 10 Hz

Q: 0.707

Appendix B | Code

Lessons

OnOpenFile

```
.....:Ventilator feedback control.bbs;.....
```

```
F = 25 ;lpm
```

```
T = 0
```

```
Array "PulseMetrics" Clear
```

```
; 0 = Pulse width
```

```
; 1 = Pulse Repetition Interval
```

```
CloseAllGraphs
```

```
Open "", r, W$
```

```
if r <> 1
```

```
    Prompt "File could not be opened!", "OK"
```

```
else
```

```
    Select Window W$
```

```
    GraphJournal New Docked
```

```
    Call "EnterFlow"
```

```
endif
```

```
End
```

```
EnterFlow
```

```
;set the flow rate
```

```
Select Window W$
```

```
RepeatFR:
```

```
GetNumber "Enter the desired flow rate: ", "lpm", F, e, false
```

```
if e <> 0
```

```
    Prompt "Exiting the program...", "OK"
```

```
    Call "ExitProgram"
```

```
else
```

```
    if F < 0
```

```
        Prompt "Need to enter Flow Rate that is > 0.", "OK"
```

```
        goto RepeatFR
```

```
    else
```

```
        Call "PulseWidth"
```

```
www.biopac.com
```

```
endif
endif
End
PulseWidth
ReEnterPulseWidth:
GetNumber "Enter Pulse Width", "sec", w, e, false
if e <> 0
    Prompt "User pressed 'Cancel or dialog error occurred. Exiting program.", "OK"
    Call "ExitProgram"
else
    if w < 0
        Prompt "Need to enter a pulse width that is > 0.", "OK"
        goto ReEnterPulseWidth
        Halt
    else
        W@ = STR$(w, 2)
        W@ = "Pulse Width entered (sec) = " + W@
        Prompt W@, "Accept", "Re-enter", "", r
        if r = 1
            Array "PulseMetrics" Set 0, w
            f = (2 / w) ;max pulse repitition frequency based on 50% duty cycle of the entered PW
            Array "PulseMetrics" Set 1, f
            Call "PulseRepFreq"
        endif
        if r = 2
            goto ReEnterPulseWidth
        endif
    endif
endif
End
PulseRepFreq
ReEnterPulseRepFreq:
Array "PulseMetrics" Get 0, p
f = 1 / (1.25 * p)
F@ = STR$(f,2)
A@ = "Enter Pulse Repitition Frequency."
A@ = A@ + "Must be less than the PRF of "
```

A@ = A@ + F@

A@ = A@ + " based on a 75% duty cycle of the entered PW."

GetNumber A@, "Hz", g, e, false

if e <> 0

 Prompt "User pressed 'Cancel or dialog error occurred. Exiting program.", "OK"

 Call "ExitProgram"

else

 if g < 0

 Prompt "Need to enter a pulse repetition frequency that is > 0.", "OK"

 goto ReEnterPulseRepFreq

 Halt

 else

 if g > f

 G@ = " " + STR\$(g,2)

 H@ = G@ + " is greater than "

 H@ = H@ + F@

 H@ = H@ + " Please re-enter a frequency that is lower."

 Prompt H@, "OK"

 goto ReEnterPulseRepFreq

 Halt

 else

 H@ = "Pulse repetition frequency entered (Hz) = " + G@

 Prompt H@, "Accept", "Re-enter", "", r

 if r = 1

 Array "PulseMetrics" Set 1, g

 Call "SetSTIM"

 endif

 if r = 2

 goto ReEnterPulseRepFreq

 endif

 endif

endif

End

SetSTIM

;relation between flowrate and voltage set to stimulator lpm --> voltage

;found through regression equation using discrete data points

Select Window W\$

$$V = (-0.0000002 * F^3) + (0.0002 * F^2) + (0.000006 * F) - 0.0115$$

Call "NewButtons"

End

StartMPAcq

Select Window W\$

Array "PulseMetrics" Get 0, G#

Array "PulseMetrics" Get 1, g

;derive the time interval between the actual pulses

$$P\# = (1 / g)$$

$$H\# = P\# - G\#$$

RemoveAllButtons

CreateButton "Done" "Done"

MP100 StartAcq

End

OnStartAcquisition

Select Window W\$

$$T\# = \text{Ticks} / 60$$

$$S\# = 0$$

MP100 Set Out0, 0

End

Done

Select Window W\$

MP100 Set Stim 0, Off

MP100 StopAcq

Call "NewButtons"

End

NewButtons

RemoveAllButtons

CreateButton "Enter Flowrate/Pulse Parameters" "EnterFlow"

CreateButton "Start Acq" "StartMPAcq"

CreateButton "Exit Program" "ExitProgram"

End

```
ExitProgram
;;AppExit
Select Window W$
CloseGraph
End
OnIdleAcquisition
if S# = 0
    if (Ticks / 60) - T# > H#
        Select Wave 8
        Edit SelectAll
        Set HCursor HCursor2 - P#, HCursor2
        X = Max
        R = X / F
        GraphJournal Append STR$(R, 2)
        GraphJournal Append "\r"
        if R < 0.1
            R = 0.1
            V = V / R
            MP100 Set Out0, V
            S# = 1
            T# = Ticks / 60
        elseif R >= 0.1
            V = V / R
            MP100 Set Out0, V
            S# = 1
            T# = Ticks / 60
        endif
    endif
elseif S# = 1
    if (Ticks / 60) - T# > G#
        MP100 Set Out0, 0
        S# = 0
        T# = Ticks / 60
    endif
endif
End
```