



CNAP[®] HD (NIBP100D-HD) IN RESEARCH



Benefits in Research
20.10.2015, V3.2

One finger sensor providing complete cardiovascular responses to tests – continuous & non invasive

SETUP & PARAMETERS





Unique combination of noninvasive continuous...

- ▶ Pulse Rate ✓
- ▶ Blood Pressure ✓
- ▶ Cardiac Output, Stroke Volume ✓
- ▶ Systemic Vascular Resistance ✓

Blood pressure, cardiac output and vascular resistance are used to differentiate between psychophysiological states.

WHY CONTINUOUS BP, CO & SVR IN RESEARCH ?

- Assessment of cardiovascular response to (psycho)physiolog. states
 - Reactions of the human body to **challenge and threat** (positive and negative stress) [1,2]



- Quantification of **emotional response** [3]

¹ Kirby, L. D., & Wright, R. A. (2003). Cardiovascular correlates of challenge and threat appraisals: A critical examination of the Biopsychosocial Analysis. *Personality and Social Psychology Bulletin*, 7, 216-233.

² Tomaka, J., Blascovich, J., Kelsey, R. M., & Leitten, C. L. (1993). Subjective, physiological, and behavioral effects of threat and challenge appraisal. *Journal of Personality and Social Psychology*, 65, 248-260.

³ Lackner, H. K., Weiss, E. M., Schuller, G., Hinghofer-Szalkay, H., Samson, A. C., & Papousek, I. (2013). I got it! Transient Cardiovascular Response to the Perception of Humor. *Biological Psychology*, 93, 33-40.

Full hemodynamics and simple setup brings huge benefits for scientific studies

FEATURES & BENEFITS OF CNAP[®] HD IN RESEARCH



SIMPLE AND QUICK SET-UP

- One finger sensor provides all parameters **noninvasively** – no placing of catheters or additional electrodes.



QUICK RECORDING

- Signals displayed only shortly after startup
- enables accurate & immediate feedback on arterial BP, cardiac output, fluid and hemodynamic status, etc.



CLINICALLY PROVEN AND VALIDATED^{1,2}

- Combination of finger sensor with NBP calibration provides high accuracy BP signal
- Continuity, accuracy & waveform dynamics are equivalent to intra-arterial measurement.
- Proven solution for consistent, repeatable results



SIMPLE DATA TRANSFER AND ANALYSIS

- Up to 4 analog output channels (BP waveform, MAP, CO*, PPV*)
- Plug & play integration into all common data acquisition systems (e.g. MP150: *AcqKnowledge*) and subject monitors.

¹ Jeleazcov et al.(2010). Precision and accuracy of a new device (CNAP) for continuous noninvasive arterial pressure monitoring: assessment during general anaesthesia. *British Journal of Anaesthesia*, 105(3):264–72.

² Wagner, J. Y. et al. Continuous noninvasive cardiac output determination using the CNAP system: evaluation of a cardiac output algorithm for the analysis of volume clamp method-derived pulse contour. *Journal of Clinical Monitoring and Computing*. doi:10.1007/s10877-015-9744-1 (2015).

*CE approval granted; FDA approval pending;

CNAP[®] HD provides easy setup for cardiac output from the same finger sensor


WAYS TO MEASURE CO IN RESEARCH

CNAP[®] HD



- Quick & simple setup
- Using a reusable finger sensor

Standard Impedance Cardiography



- Complex and time-consuming setup requires to undress upper body
- Placing of disposable electrodes

CNAP[®] HD combines three core technologies

HOW DOES CNAP[®] HD WORK?



COMBINATION OF 3 TECHNOLOGIES/METHODS:

- 1. Vascular unloading principle**
 - Blood pressure, pulse rate from finger sensor
- 2. Standard upper arm cuff for calibration (NBP)**
 - Reference BP value at heart level for automatic scaling
- 3. Pulse Contour Analysis (PCA)**
 - SV, CO, SVR





Integrated pressure chambers measure blood pressure continuously



CONTINUOUS BLOOD PRESSURE

- 1 Infrared light sensors pick up blood volume and flow.
- 2 In- and deflating pressure chambers follow blood pulsation and keep blood flow constant. The resulting pressure in the finger sensor corresponds to the real arterial pressure.¹
- 3 The pressure system is controlled by multiple digital feedback loops and the “VERIFI-algorithm” for high fidelity signal processing and artifact rejection.²




¹ Jan Penaz, Brno, Czech Republic 1973

² Fortin J et al.: Continuous non-invasive blood pressure monitoring using concentrically interlocking control loops. Comput in Biol Med 36 (2006) 941 – 957

Interaction of fingersensor with NBP provides high accuracy BP signal – equivalent to invasive arterial line

CALIBRATION TO CLINICAL GOLD STANDARD




Cardiac output measure obtained from accurate blood pressure waveform using well-known pulse contour analysis

CNAP[®] HEMODYNAMICS

➤ Applying Ohm's Law to Hemodynamics


Voltage = Current x Resistance

$$\text{MAP} - \text{CVP} = \text{CO} \times \text{SVR}$$



$$\text{MAP} \approx \text{CO} \times \text{SVR}$$

MAP ... mean arterial pressure
CVP ... central venous pressure*
CO ... cardiac output
SVR ... systemic vascular resistance



▶	MAP	derived from accurate blood pressure waveform provided by Vascular unloading
▶	CO	is calculated from pulse contour analysis (PCA)
▶	SVR	results from CO and MAP
▶	SVI, CI, SVRI	absolute values can be indexed to body surface area


* CVP ≈ 7 [mmHg] and can be neglected

CNAP[®] arterial pressure is comparable with an invasive arterial line, even during general anesthesia

CLINICAL VALIDATION: CNAP[®] BLOOD PRESSURE

Erlangen, Germany^[1]:

- 88 Patients (Neuro, Abdominal, Cardiac Surgery)
- Bland-Altman Plot for MAP:
 - Bias (SD) = -1.6 (11.0) mmHg
- Precision (trending):
 - CNAP[®] as accurate as invasive BP (3,2 mmHg for MAP)
- 82.1% fast blood pressure changes *detected simultaneously* by CNAP[®] and IBP





¹ Jeleazcov et al. Precision and accuracy of a new device (CNAP) for continuous non-invasive arterial blood pressure monitoring: assessment during general anesthesia. British Journal of Anesthesia 2010 vol. 105 (3) pp. 264-72

CO with CNAP[®] HD is inter-changeable with clinical standards and tracks changes reliably

CLINICAL VALIDATION: CARDIAC OUTPUT WITH CNAP[®] HD

- **HIGH ACCURACY** of non-invasive CO with CNAP[®] HD compared to invasive transpulmonary thermodilution¹. Bland-Altman plots show:
 - Small bias (0.2 l/min) and narrow limits of agreement (LoA: -1.7 to 2.2)
 - Percentage error: 25% (according to strict Critchley standard of $\leq 30\%$ ^{2,3})
- **HIGH TRENDING ABILITY** of non-invasive CO with CNAP[®] HD compared to invasive clinical standards
 - Concordance rates $\geq 90\%$ ⁴ (according to standard⁵)

CO with CNAP[®] HD vs. Thermodilution



¹ Wagner, J. Y., Grond, J., Fortin, J., Negulescu, I., Schöfthaler, M., & Saugel, B. (2015). Continuous noninvasive cardiac output determination using the CNAP system: evaluation of a cardiac output algorithm for the analysis of volume clamp method-derived pulse contour. *Journal of Clinical Monitoring and Computing*. doi:10.1007/s10877-015-9744-1

² Critchley, L. A., & Critchley, J. A. (1999). A meta-analysis of studies using bias and precision statistics to compare cardiac output measurement techniques. *Journal of clinical monitoring and computing*, 15(2), 85–91

³ Peyton PJ, Chong SW: Minimally invasive measurement of cardiac output during surgery and critical care. *Anesthesiology* 2010; 113:1220-35.

⁴ CNAP[®] HD validation data for CE mark and prepared for publication in peer-reviewed journal – data available upon request

⁵ Critchley LA, Lee A, Ho AMH: A critical review of the ability of continuous cardiac output monitors to measure trends in cardiac output. *Anesth Analg* 2010; 111:1180-92.

CNAP[®] is used in a variety of applications & scientific studies

RESEARCH EXAMPLES WITH CNAP[®]

BIOPAC Customers

Psychophysiology: Challenge & Threat

Cornick, Blascovich et al. from the University of California in Santa Barbara investigated the cardiovascular consequences of exercise in obese subjects in a stressful environment. Cardiovascular and self-report measures were recorded. Findings indicated that those who are more self-aware showed cardiovascular response patterns indicative of threat.¹

Psychophysiology: Emotional Response

At Pitzer College, California, Noble et al. performed a study to assess the effects of mainstream media on women's physiological and psychological functioning.²

The Naval Aerospace Medical Research Laboratory studied the **detection of deception** by use of continuous blood pressure³.

Sanchez-Gonzalez et al. evaluated the cardiovascular reactivity to anxiety in young and middle-aged individuals⁴.

Psychophysiology: Interaction Hemodynamics & Brain

The Institute for Knowledge Discovery from Graz, Austria studied the influence of slow blood pressure oscillations on self-paced movements. "Free will" is not only brain directed but also a matter of blood pressure.⁵

¹ Cornick, J. E., & Blascovich, J. (2015). Consequences of objective self-awareness during exercise. *Health Psychology Open*, 2(2), 2055102915598088. doi:10.1177/2055102915598088

² Noble, M. L. (2012). The Effect of Mainstream Media on Body Image and Stress Reactivity in Latina Females. *Pitzer SeniorTheses*.

³ Taylor, M. K., Horning, D. S., Chandler, J. F., Phillips, J. B., Khosravi, J. Y., Bennett, J. E., ... Gao, H. (2011). A Comparison of Approaches To Detect Deception. *Technical Report, Naval Aerospace Medical Research Laboratory*, (11).

⁴ Sanchez-Gonzalez, M. A., Guzik, P., May, R. W., Koutnik, A. P., Hughes, R., Muniz, S., ... Fincham, F. D. (2015). Trait anxiety mimics age-related cardiovascular autonomic modulation in young adults. *Journal of Human Hypertension*, 29(4), 274–80. doi:10.1038/jhh.2014.72

⁵ Pfurtscheller, G., Ortner, R., Bauernfeind, G., Linortner, P., & Neuper, C. (2010). Does conscious intention to perform a motor act depend on slow cardiovascular rhythms? *Neuroscience Letters*, 468(1), 46–50

CNAP[®] is used in a variety of applications & scientific studies

RESEARCH EXAMPLES WITH CNAP[®]

BIOPAC Customers

Physiology

Hurr et al. from the Department of Kinesiology and Health Education (The University of Texas at Austin; USA) determined whether there is a difference in cerebral vasodilatory capacity in response to rebreathing-induced hypercapnia between African Americans and Caucasian Americans.⁶

Sports

Telles et al. assessed heart rate variability, non-invasive arterial BP, and respiration rate, during alternate nostril yoga breathing and breath awareness sessions.⁷

Gerontology

Pereira et al. from Marquette University, USA, investigated the steadiness of muscle contractions in young and old adults during high and low cognitive demand tasks.⁸

Contact us for a complete bibliography of CNAP[®] / NIBP100D studies

⁶ Hurr, C., Kim, K., Harrison, M. L., & Brothers, R. M. (2015). Attenuated cerebral vasodilatory capacity in response to hypercapnia in college-aged African Americans. *Experimental Physiology*, 100(1), 35–43. doi:10.1113/expphysiol.2014.082362

⁷ Telles, S., Sharma, S. K., & Balkrishna, A. (2014). Blood Pressure and Heart Rate Variability during Yoga-Based Alternate Nostril Breathing Practice and Breath Awareness. *Medical Science Monitor Basic Research*, 20, 184–93. doi:10.12659/MSMBR.892063

⁸ Pereira, H. M., Spears, V. C., Schlinder-Delap, B., Yoon, T., Nielson, K. A., & Hunter, S. K. (2015). Age and sex differences in steadiness of elbow flexor muscles with imposed cognitive demand. *European Journal of Applied Physiology*, 115(6), 1367–79. doi:10.1007/s00421-015-3113-0