


# BIO PAC




**BIOPAC**  
Systems, Inc.  
Registered to ISO 9001:2008

*Inspiring people and  
enabling discovery about life*

# Housekeeping

- Attendees are on Mute
- Headset is Recommended!
- Questions addressed at end of Webinar and in a Q&A follow-up document



## Your Participation

Open and hide your control panel

Join audio:

- Choose “Mic & Speakers” to use VoIP
- Choose “Telephone” and dial using the information provided

Submit questions and comments via the Questions panel

**Note:** Today’s presentation is being recorded and will be provided when available.



# BIOPAC Guide to EEG II



**Frazer Findlay**


CEO, BIOPAC Systems, Inc.



**Chris Berka**

CEO, Advanced Brain Monitoring





Audience Poll





Wireless, Wearable Data Recording

# PHYSIOLOGY ANYWHERE

Complete Solutions for Life Science Data Acquisition & Analysis.

LEARN MORE







Over 97% of top universities run BIOPAC Systems

# THE WORLD DISCOVERS WITH BIOPAC

Solutions for life science research and education

# THE WORLD DISCOVERS WITH BIOPAC


**BIOPAC systems are cited  
over 27,000**





# THE WORLD DISCOVERS WITH BIOPAC

BIOPAC Systems cited ~4,000 times  
for EEG





# INTRODUCTION

Wireless EEG with B-ALERT

Subject Setup

Recording wireless EEG

Recording wireless EEG with other signals


Cognitive States Metrics - background

Recording Cognitive States Metrics and Analysis

Summary

Q&A





# Demonstration of *AcqKnowledge* Software



# Combining EEG with other Physiological Signals



## **MP160 System**

16-Channels of  
Tethered or  
Wireless Signals



## **B-Alert System**

9-Channels of EEG  
and Cognitive  
States Metrics

# B-Alert Wireless EEG System

## B-ALERT<sup>®</sup>X10




## ACQKNOWLEDGE<sup>™</sup>

### **B-ALERT X10**

- Nine channel EEG
- Wireless - Bluetooth
- Combine with MP160 for additional channels of data
- Uses linked mastoids for reference signal
- Built-in accelerometer

# B-Alert Wireless EEG System




## B-ALERT X10

- Locations F3, FZ, F4, C3, CZ, C4, P3, POz, P4
- Leads for ECG
- Synchronize with MP160 using CBLX10 Cable



# Correct Fit

Nasion to Occipital



Crus of Helix



**B-Alert Sensor Strip Sizing Chart**



S = Small M = Medium

Nasion to Occipital		Crest of Helix to Crest of Helix																			
		24.5		25.5	26	26.5	27	27.5	28	28.5	29	29.5	30	30.5		31.5	32	32.5	33	33.5	34
	30.5	XS	XS	XS	XS	XS	XS	XS	XS	XS	XS	XS	-	-	-	-	-	-	-	-	-
	31	XS	XS	XS	XS	XS	XS	XS	XS	XS	XS	XS	-	-	-	-	-	-	-	-	-
	31.5	XS	XS	XS	XS	XS	XS	XS	XS	XS	S	S	S	S	S	S	S	S	-	-	-
	32	XS	XS	XS	XS	XS	XS	XS	XS	XS	S	S	S	S	S	S	S	S	-	-	-
	32.5	XS	XS	XS	XS	XS	XS	S	S	S	S	S	S	S	S	S	S	S	M	M	M
	33	XS	XS	XS	XS	XS	S	S	S	S	S	S	S	S	S	M	M	M	M	M	M
	33.5	XS	XS	XS	XS	S	S	S	S	S	S	S	S	S	S	M	M	M	M	M	M
	34	XS	XS	XS	S	S	S	S	S	S	S	S	S	S	S	M	M	M	M	M	M
	34.5	XS	XS	XS	S	S	S	S	S	S	S	S	M	M	M	M	M	M	M	M	M
	35	XS	XS	XS	S	S	S	S	S	S	S	S	M	M	M	M	M	M	M	M	M
	35.5	XS	XS	XS	S	S	S	S	S	S	S	S	M	M	M	M	M	M	M	M	M
	36	XS	XS	S	S	S	S	S	S	S	S	M	M	M	M	M	M	M	M	M	M
	36.5	-	S	S	S	S	S	S	S	S	S	M	M	M	M	M	M	M	M	M	M
	37	-	-	-	S	S	S	S	S	S	M	M	M	M	M	M	M	M	M	M	M
37.5	-	-	-	S	S	S	S	S	S	M	M	M	M	M	M	M	M	M	M	M	

XS = extra small strip available on request

- = head dimensions outside the suitable range for ABM strips

# B-Alert Setup



## B-ALERT X10

- Electrode array
- Adhesive foam pads
- Fill cavity with gel
- Apply array to the subject




Fuzzy side of neoprene



Triangle Strip Connector




Mastoid Leads



Audience Poll




# B-Alert Setup



## B-ALERT X10

- Electrode array
- Adhesive foam pads
- Fill cavity with gel
- Apply array to the subject

## Subject Preparation



Fuzzy side of neoprene





Triangle Strip Connector



Mastoid Leads

# Impedance Check





# Demonstration of *AcqKnowledge* Software with B-Alert

# Combining EEG with other Physiological Signals



**MP160 System**

16-Channels of  
Tethered or  
Wireless Signals



**CBLX10**

Isolated  
Synchronization  
Cable




**B-Alert System**

9-Channels of EEG  
and Cognitive  
States Metrics






# Demonstration of *AcqKnowledge* Software with B-Alert




Audience Poll



## The B-Alert Cognitive State Metrics<sup>1</sup>

Four (4) EEG-based classifications that span a continuum from highly engaged and processing sensory inputs, to fatigued and inattentive.




<sup>1</sup> Johnson, R. R., D. P. Popovic, et al. (2011). Drowsiness/alertness algorithm development and validation using synchronized EEG and cognitive performance to individualize a generalized model. *Biol Psychol* 87(2): 241-250 .




# The Cognitive Workload Metric<sup>2</sup>

In training environments, the most effective programs sustain optimal workload for learning, and avoid extreme measures.



<sup>2</sup> Berka C, Levendowski DJ, Lumicao MN, Yau A, Davis G, et al (2007). EEG correlates of task engagement and mental workload in vigilance, learning, and memory tasks. *Aviat. Space Environ. Med.* 78, B231–B244. May 2007.

# B-Alert Cognitive States Metrics



## B-ALERT X10

### Metrics

- Fatigue
- Stress
- Confusion
- Engagement
- Workload


### Regions

- Executive functions (prefrontal cortex)
- Motor control (motor-strip)
- Visual system activity (parietal region)

The background image shows the AcqKnowledge software interface. At the top, there is a toolbar with various icons for file operations, settings, and data manipulation. Below the toolbar, a message bar displays "Message: Image 8 Puma". The main area contains several channels of data plotted on a grid. The top channel is a red waveform showing high-frequency, regular oscillations. The middle channel is a blue waveform showing a low-frequency signal with a prominent peak. The bottom channel is a purple waveform showing a low-frequency signal with some sharp spikes. The text "Demonstration of AcqKnowledge Software with B-Alert" is overlaid on the bottom half of the image in white font.

# Demonstration of *AcqKnowledge* Software with B-Alert





Stimulus Presentation with  
Triggers for Synchronization

SuperLab® - E-Prime®  
Presentation®

Virtual Reality

BIOPAC VR Solutions  
offer real-time feedback  
loops and event marking

# Stimulus Presentation Options


STP100C



ESU







# SUMMARY

9 Channel wireless EEG

Sizing the subject for the correct array

Linked Mastoids for reference

Cognitive States Metrics

Subject specific Definitions file

Synchronization with other signals



# Q&A Session:

Please submit questions for our guest speakers through the Questions Window. While all questions cannot be answered during our live session, all will be reviewed and answered following our event.



...Ask about the T4 conference



# BIOPAC Physiology T4 Conference



- Meet us in Santa Barbara!
- July 24-26, 2017
- University of California at Santa Barbara, USA
- Over 50 Human Physiology workshops!
- World-renown researchers present
- [www.biopac.com](http://www.biopac.com)



**For more information:**

[www.biopac.com](http://www.biopac.com)  
[info@biopac.com](mailto:info@biopac.com)

View upcoming and on-demand webinars at:  
[www.biopac.com/webinars](http://www.biopac.com/webinars)

**Thank you for your time and attention!**

Please complete our survey upon  
exiting the webinar.  
Thank you for attending!