

MAXACQ-M *Discontinued July 2015*

This software interface will permit users to work with MP Research Systems for Mac (MP150WS or MP36RWS) to collect data into Max/MSP/Jitter via Network Data Transfer protocol; compatible with Mac OS X 10.5-10.6.

MaxAcq_M is a Java (mxj) object for Max/MSP/Jitter for communicating with BIOPAC's *AcqKnowledge* server. It interfaces with *AcqKnowledge* server's streaming capability, outputting data channels as Jitter matrices.



Prerequisites: *Discontinued July 2015—not validated for OS or MAX updates*

- *AcqKnowledge* software with network capability from BIOPAC
 - Max5 from Cycling '74 (www.cycling74.com)
 - Java runtime environment

Workflow:

1. Open *AcqKnowledge* application.
2. Open Max5 and a patcher with the MaxAcq mxj object. MaxAcq will automatically detect the *AcqKnowledge* server if it exists.
3. Enable sensor channels in *AcqKnowledge* either by loading in a template file through the MaxAcq 'loadTemplate' message or directly within the *AcqKnowledge* application
4. Enable channels for streaming (e.g. calc 0, analog 0, digital 0) with the MaxAcq 'stream' message.
5. Send a stream of bang messages to MaxAcq (usually with a metro) to output data as it arrives.

Optional Support License

BIOPAC Developer products are intended for users with some programming knowledge. Issues that the Support Department can resolve do not require and are not counted toward a license. Optional Support Licenses provide for up to 5, 10, or 20 software-related programming issues that require a programmer to resolve. Click for [Support License](#) options.

In addition to the Help file, a bare-bones Max example for testing purposes is included in examples/example.test.maxpat.

Test Run:

1. Open *AcqKnowledge* application; the program must be running for MaxAcq to detect it as a server.
2. Open Max5 and within Max5, open the file: example.test.maxpat, found at Applications/Max5/examples/MaxAcq/example.
3. Within the "example.test.maxpat" click on the button above the 'opendialog' text box to load an *AcqKnowledge* template file located at: /Applications/Max5/examples/MaxAcq/template/acq_template.acq.
4. Within the "example.test.maxpat" window, under text "Output Data," click on the toggle below it to enable output to the patch window (example.test.maxpat file window).
5. Click on the 'stream calc 0' message. The *AcqKnowledge* graph file should start acquiring data and the amplitude data from calculation channel CH0 should be streamed to and printed within the patch window.