

MOBITA WEARABLE, WIRELESS PHYSIOLOGICAL MONITORING & LOGGING

- 32 channels of data in the palm of your hand
- Mobile physiologic measurement with Mobita & AcqKnowledge

MOBITA-W SERIES

Mobile 32-Channel Biopotential systems



MOBITA-W-32EEG System includes ConfiCap MB-32EEG-CAP-A

MOBITA-W-20EEG System includes ConfiCap MB-20EEG-CAP-B

MOBITA-W-12+20 System includes ConfiCap MB-12+20-CAP

MB-CAP SERIES

Mobita ConfiCap only



MB-32EEG-CAP-A

MB-20EEG-CAP-B

MB-12+20-CAP

Mobita Overview

Mobita® is a new wearable physiological signal amplifier system that can record up to 32 channels of high-fidelity wireless biopotential data, including ECG, EEG, EGG, EMG, and EOG data. The device includes a trigger channel that can be used to synchronize the system with other devices or data streams. When the onboard accelerometer is used with AcqKnowledge's Actigraphy feature, it is possible to evaluate a subject's activity levels.

Record in the lab or out, from simple to tough and demanding measurement situations. The system can either telemeter data back to a computer running AcqKnowledge for real-time display and analysis of the signals, or record it locally for later download. Easily switch between live or logging modes to suit your research protocol. Compact design, powerful specs, and impressive flexibility combine to create the ultimate solution for mobile physiologic measurement!

Key System Features

- 32 channels of wireless biopotential data
- Fully-integrated in AcqKnowledge® software
- Built-in 3D accelerometer for position information
- Trigger channel supports TTL inputs from 3rd - party hardware
- Built-in WiFi telemetry (range typical > 10 m indoors)
- No cable movement artifacts
- No filtering (including Notch filter) for true unadulterated signal quality
- True DC recording
- 24 bit data resolution
- Flash disk recording (up to 16 GB) for data back-up and holter applications
- Rechargeable Li-Po battery
- Rugged construction: sturdy, dustproof enclosure

Featured Applications

- Psychophysiology
- Neuroscience
- Exercise Physiology & Sports
- Gait & Movement Analysis
- Brain Computer Interface
- Sleep Studies
- Neuromarketing
- Home-based Ambulatory Testing

Flexibility

With AcqKnowledge and Mobita®, the system is quickly configured to do the work of multiple systems without the added cost of multiple amplifiers.

Quickly change the electrode configuration or signal type by swapping out the ConfiCap. Simply disconnect one header and snap on a new configuration for a completely different application. Record a 12-lead ECG while recording EEG and EMG data all with the same device.

Analyze with AcqKnowledge

- Powerful automated analysis routines for ECG, HRV, EEG, EMG, EGG, and many more!
- Intuitive user interface with fully customizable display

- Video Tutorials on key features and analysis routines
- Guided channel and acquisition setup with presets and quickstarts

MOBITA-W System Options



MOBITA-W-32EEG System
 32-Channel EEG



MOBITA-W-20 EEG System
 10/20 EEG + 13 open leads



MOBITA-W-10+12 System
 32-Channel Biopotentials

Complete systems include one Mobita hardware unit, one ConfiCap configuration with accessories, and AcqKnowledge software. Each flexible Mobita system records up to 32 channels of data at up to 2K s/s. Stream data live into AcqKnowledge or log data for later upload. Swap ConfiCaps to change experiment protocols.

- **MOBITA-W-32EEG System** This system includes one Mobita hardware unit and one ConfiCap MB-32EEG-CAP-A with medium 32-channel headcap (H2O-CAP-M; other cap sizes available) and water electrodes. Water electrodes eliminate the need for gel. Individually shielded inputs.
- **MOBITA-W-20EEG System** This system includes one Mobita hardware unit and one ConfiCap MB-20EEG-CAP-B configured for 10/20 EEG; and EEG Cap (medium CAP100C; other sizes available). Leaves 13 open TouchProof inputs for user’s choice of biopotential data. Unshielded inputs.
- **MOBITA-W-12+20 System** This system includes one Mobita hardware unit and one ConfiCap MB-12+20-CAP to record up to 32 channels of biopotential data. Unshielded inputs.

Specifications

Sampling

Resolution: 24.414 nV/bit, referred to input
 Sampling rate: 2000, 1000, 500, 250 Hz
 Channel bandwidth: DC up to 0.2 x sample freq

Input:

Input signal difference: 409.6 mV pp
 Input common mode range: -2.0 V . +2.0 V
 Gain factor: 10
 Noise: < 0.4 µV RMS @ 0.1 . 10 Hz
 Input Impedance: >10 G
 CMRR: > 100 dB typical
 # of channels: Up to 32 analog
 Power supply: Battery Li-Polymer with protection circuit

Accelerometer

Range: ± 16 g
 Resolution: 13 bit

Sensitivity: 3.9 mg/bit

Battery life:

WiFi mode: 8-10 hours (environment dependent)
 Logger Mode: 17-19 hours
 Hybrid Mode: WiFi and Logger: 8-10 hrs (environment dependent)

Note: The number of channels enabled does NOT significantly influence battery life.

Filtering

Filter: No filtering within channel bandwidth
 Type: Unipolar, bipolar (user configurable from unipolar inputs)
 Trigger: Either trigger (TTL) or generic (e.g., RS232 compatible) digital inputs possible through custom designed ConfiCap

Dimensions: 150 x 70 x 25 mm (with ConfiCap)

Computer Requirements

Computer should be running Windows 7 64-bit or Windows 8 64-bit with a Core i5 or a Core i7 processor.


Notes ! No support provided for operating systems older than Windows 7.

Slower computers may be able to use WiFi mode with the Mobita, but it may not be possible to transfer or import the logged data.

ConfiCap™ Options


ConfiCaps allow you to quickly change the configurations of the inputs for specific applications (i.e., 32-ch EEG, EMG, or combinations of ECG/EMG/EEG, etc.), customize your own design, or create protocol driven configurations when recording with a MOBITA-W System. Each channel is unipolar (single-ended) and AcqKnowledge is easily configured to create unique montages.

- **MB-32EEG-CAP-A 32-Channel EEG ConfiCap**



The MB-32EEG-CAP-A is a complete assembly for the Mobita wearable biopotential system that interfaces with a 32-channel electrode cap. This particular EEG cap uses water electrodes, which eliminates the need for gel. The assembly also includes a trigger connector for synchronization with other devices. The electrodes terminate in a Mobita confi-cap connector. Snap the assembly to the Mobita unit and attach the cap to a subject to record 32 channels of data for either in laboratory telemetry or remote data logging applications.

TIP: Use a permanent marker on the white top of each input to indicate its montage position to help expedite setup.



H2O-CAP Series Headcaps for Mobita Water Electrodes

- H2O-CAP-SMALL** 50-54 cm
- H2O-CAP-MEDIUM** 54-58 cm
- H2O-CAP-LARGE** 58-62 cm

Headcaps include 32 grommets for Mobita water-based electrodes. Each MB-32EEG-CAP-A assembly includes one user-specified H2O-CAP size); individual headcaps can be used to add to or replace the included cap.

- **RX-H2O-ELECT Replacement Water Electrode Pack**

Replacement pack of 100 water electrode strips for the Mobita EEG system. These small strips of absorbent material must be rolled up and inserted into the electrode cavity and placed on the EEG cap prior to recording. The Mobita wearable physiological signal amplifier system records 32 channels of high-fidelity wireless EEG data and includes 100 water electrode inserts. This pack of 100 strips can be used for additional subjects.



• **MB-20EEG-CAP-B EEG 10/20 + 13 TP Adapters**



The MB-20EEG-CAP-B is a complete assembly for the Mobita wearable biopotential system that interfaces with a 10/20 electrode cap and 1.5 mm Touchproof sockets for adding additional signals.

This combination interface allows for a full 10/20 EEG, plus optional biopotential signals for EOG, EMG, and ECG. Snap the assembly to the Mobita unit and attach the cap to a subject to record 32 channels of data for either in laboratory telemetry or remote data logging applications.

• **MB-12+20-CAP 12 Surface Electrodes + 20 TP Adapters**



The MB-12+20-CAP is a complete assembly for the Mobita wearable biopotential system that interfaces with 12 snap fit electrode leads and 20 Touchproof (1.5 mm) sockets. Record 32 channels of biopotential data using a variety of electrode configurations including both disposable and reusable options. Connect to the Mobita and the subject and record up to 32 channels of data for either in laboratory telemetry or remote data logging applications.

Breakout boxes are available. Contact [BIOPAC](http://www.biopac.com) for more information.