


Application Note 271 Synchronize B-Alert and MP150 systems running on separate computers


1. Connect the hardware as shown in the figure below.



Attach to the EKG (ECG) input of the B-Alert System

Attach the 3.5 mm phone jack connector to the Analog Out (A0) on the front of the UIM100C

2. Set up the **MP150 System** and *AcqKnowledge* software as required for the experiment.
3. In the *AcqKnowledge* software for the **MP150 System**:
 - Choose **MP150 > Set Up Stimulator** and set parameters as shown in the following figure:



4. In the *AcqKnowledge* software for the **B-Alert System**:
 - Choose **B-Alert > Set Up Channels**.
 - Check **Acquire and Plot** for Analog Channel 1 (A1/EKG).
5. Click the **Start** button in the *AcqKnowledge* Software for the **B-Alert System**.

Input channels setup for 'B-Alert 300103014/10ch/16bits'




Acquire	Plot	Value	Channel	Label
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A1	EKG
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A2	Poz
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A3	Fz
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A4	Cz

IMPORTANT! You must start *AcqKnowledge* for the B-Alert System before starting *AcqKnowledge* for the MP150 system.

6. Click the **Start** button in the AcqKnowledge software for the **MP150 System**.
7. Observe the **B-Alert** data. A negative and positive spike* should appear and resemble the figure below. (It may be helpful to autoscale CH1 and hide additional channels if necessary for a better view.)
 - The event marker shows the beginning of the stimulus. (This marker was added manually for reference.)

* Due to filtering applied by the EKG channel, the square wave will appear truncated in the B-Alert waveform. If desired, an AC-coupled sine wave can be used instead for a more accurate representation.

8. After data recording is complete, return to the beginning of the **B-Alert** graph where the stimulus occurred.
9. Highlight all data before the stimulus, and choose **Edit > Clear**.



10. **MP150** Data is now synced to **B-Alert** Data (time point 'zero' of the **MP150** data file equals time point 'zero' of the **B-Alert** data file).

