

## V. DATA ANALYSIS

### FAST TRACK Data Analysis

1. Enter the **Review Saved Data** mode.

- Note Channel Number (CH) designations:

*Channel*     *Displays*

**CH 1**     **EEG**

**CH 40**     **alpha**

**CH 41**     **alpha RMS**

- Note measurement box settings:

*Channel*     *Measurement*

**CH 1**     **Stddev**

**CH 40**     **Stddev**

**CH 41**     **Mean**


**CH 40**     **Freq**

2. Set up your display window for optimal viewing of the entire recording.

### Detailed Explanation of Data Analysis Steps

If entering **Review Saved Data** mode from the Startup dialog or lessons menu, make sure to choose the correct file.

The data should resemble Fig. 4.13.



**Fig. 4.13 Example data**

The measurement boxes are above the marker region in the data window. Each measurement has three sections: channel number, measurement type, and result. The first two sections are pull-down menus that are activated when you click them.


#### Brief definition of measurements:

**Stddev:** Standard deviation is a measure of the variability of data points. The advantage of the Stddev measurement is that extreme values or artifacts do not unduly influence the measurement.

**Mean:** Displays the average value in the selected area.

**Freq:** Converts the time segment of the selected area to frequency in cycles per second

The selected area is the area selected by the **I-beam** tool (including endpoints).

**Note:** The append event markers  mark the beginning of each recording. Click on (activate) the event marker to display its label.

#### Useful tools for changing view:

**Display menu:** Autoscale Horizontal, Autoscale Waveforms, Zoom Back, Zoom Forward

**Scroll Bars:** Time (Horizontal); Amplitude (Vertical)

**Cursor Tools:** Zoom Tool

**Buttons:** Overlap, Split, Show Grid, Hide Grid, -, +

**Hide/Show Channel:**  $\delta$ Alt + click (Windows) or  $\delta$ Option + click (Mac) the channel number box to toggle channel display.

**Data Analysis continues...**

- Use the I-Beam cursor to select the first data recording.



A

- Repeat the measurements for each of the data recordings.



A


- Zoom in on a small section of the Recording 1 data.
- Use the I-Beam cursor to select an area from one peak to the next in the **alpha** band (CH 40).



B

- Answer the questions at the end of the Data Report.
- Save** or **Print** the Data Report.
- Quit** the program.


**END OF DATA ANALYSIS**



**Fig. 4.14 First data recording selected**

Be sure to zoom in far enough so that you can easily measure the frequency of the **alpha** wave.

Fig. 4.15 shows a sample setup for measuring the frequency in the **alpha** band (CH 40).



**Fig. 4.15 Alpha wave frequency measurement**

An electronically editable **Data Report** is located in the journal (following the lesson summary,) or immediately following this Data Analysis section. Your instructor will recommend the preferred format for your lab.

### END OF LESSON 4

Complete the Lesson 4 Data Report that follows.

# ELECTROENCEPHALOGRAPHY II

• EEG II

## DATA REPORT

Student's Name: \_\_\_\_\_

Lab Section: \_\_\_\_\_

Date: \_\_\_\_\_

### I. Data and Calculations

#### Subject Profile

Name: \_\_\_\_\_

Height: \_\_\_\_\_

Age: \_\_\_\_\_

Gender: Male / Female

Weight: \_\_\_\_\_

#### Amplitudes

A. Complete Table 4.1 with the amplitudes of the recorded data in the control and experimental conditions. Calculate the difference for the Alpha-RMS Mean between the Experimental Conditions and the Control, and then summarize whether the Experimental Mean was larger (+), smaller (-), or the same (=) as the Control Mean.

**For example:** To calculate Alpha-RMS Difference for the "Mental Arithmetic" recording, subtract the "Eyes Closed (Control)" Alpha-RMS value from the measured "Mental Arithmetic" Alpha-RMS value.

Table 4.1

Condition	EEG		Alpha		Alpha-RMS		Alpha-RMS Difference (Exp. - Control)	Alpha-RMS Summary (+, -, =)
	1	Stddev	40	Stddev	41	Mean		
Eyes closed (Control)								
Mental arithmetic								
Recovering from hyperventilation								
Eyes open								

#### Frequency

B. What is the frequency of an alpha rhythm from "Eyes closed" data? 40 Freq = \_\_\_\_\_ Hz

Does this agree with the expected values? Yes No

### II. Questions

C. Refer to Table 4.1: When was the general amplitude of the EEG highest?

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

D. Refer to Table 4.1: When were the alpha wave levels highest?

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

E. Refer to Table 4.1: How do your results compare with the information presented in the Introduction?

---

---

F. Did Subject need to concentrate during math problems? Yes No  
How would the level of concentration required affect the data?

---

---

---

---

G. What might account for the amplitude difference of waves recorded from a subject tested alone, in a darkened room, and subjects tested in a lab full of students?

---

---

---

H. Which conditions produced the lowest alpha activity?

---

---

---

**III. OPTIONAL Active Learning Portion**

A. *Hypothesis*

---

---

---

---

B. *Materials*

---

---

---

---

---

---

C. *Method*

---

---

---

---

---

---

D. *Set Up*

---

---

---

---

---

---

E. *Experimental Results*

---

---

---

---

---

---