

V. DATA ANALYSIS

FAST TRACK Data Analysis

- Enter the **Review Saved Data** mode.

- Note channel number (CH) designations:

<i>Channel</i>	<i>Displays</i>
CH 1	Airflow (hidden)
CH 2	Volume

- Note the measurement box settings:

<i>Channel</i>	<i>Measurement</i>
CH 2	P-P
CH 2	Max
CH 2	Min
CH 2	Delta

Data Analysis continues...

Detailed Explanation of Data Analysis Steps

If entering **Review Saved Data** mode from the Startup dialog or Lessons menu, make sure to choose the correct file.




Fig. 12.23 Example data

All measurements will be performed on the Volume (CH 2) data. The Airflow (CH 1) data, used to calculate volume, is hidden to avoid confusion. It can be shown by ⌘Alt + click (Windows) or ⌘Option + click (Mac) the channel number box.

The measurement boxes are above the marker region in the data window. Each measurement has three sections: channel number, measurement type, and result. The first two sections are pull-down menus that are activated when you click them.

Brief definition of measurements:

P-P (Peak-to-Peak): Subtracts the minimum value from the maximum value found in the selected area.

Max: Displays the maximum value in the selected area.

Min: Displays the minimum value in the selected area.

Delta: Computes the difference in amplitude between the last point and the first point of the selected area.

The selected area is the area selected by the I-Beam tool (including endpoints).

Useful tools for changing view:

Display menu: Autoscale Horizontal, Autoscale Waveforms, Zoom Back, Zoom Forward

Scroll Bars: Time (Horizontal); Amplitude (Vertical)

Cursor Tools: Zoom Tool

Buttons: Overlap, Split, Show Grid, Hide Grid, -, +

Hide/Show Channel: ⌘Alt + click (Windows) or ⌘Option + click (Mac) the channel number box to toggle channel display.

2. Review the measurements described in the Introduction to identify the appropriate selected area for each.




Fig. 12.24 Measurement areas for respiratory volumes and capacities

3. Calculate the Predicted Vital Capacity, then measure VC and then compare the two.



The selected area should start just prior to the maximum peak and end just after the minimum peak. The P-P (peak to peak) measurement displays the VC.




Fig. 12.25 Example selected area; P-P measures VC

4. Take two measures on the third TV cycle:

a) Use the **I-beam** cursor to select the **inhalation** of cycle 3 and note the P-P result (Fig. 12.26). The selected area should be from the valley to the peak of the third cycle.



The P-P measurement in Fig. 12.26 represents the first value required for the averaged TV calculation.




Fig. 12.26 Example of cycle 3 – Inhale selection to measure TV

b) Use the **I-beam** cursor to select the **exhalation** of cycle 3 and note the P-P result (Fig. 12.27). The selected area should be from the peak to the valley of the third cycle.



The P-P measurement in Fig. 12.27 represents the second value required for the averaged TV calculation.




Fig. 12.27 Example of cycle 3 – Exhale selection to measure TV

Data Analysis continues...

- Repeat TV measurements, as in Step 4, but on cycle 4 data. Calculate average value of all four TV measurements.



B

- Use the I-beam cursor and measurement tools to record the volumes and capacities required by the data report (defined in Fig. 12.24).



B

- Answer the questions at the end of the Data Report.
- Save** or **Print** the data file.
- Quit** the program.

Note that the Delta measurement requires precise placement of the selected area.




Fig. 12.28 Example selection for measurements of TLC (Max) and IRV (Delta)

An electronically editable **Data Report** is located in the journal (following the lesson summary,) or immediately following this Data Analysis section. Your instructor will recommend the preferred format for your lab.

END OF DATA ANALYSIS

END OF LESSON 12

Complete the Lesson 12 Data Report that follows.

PULMONARY FUNCTION I

- *Volumes and Capacities*

DATA REPORT

Student's Name: _____

Lab Section: _____

Date: _____

Subject Profile

Name: _____ Height: _____ Gender: Male / Female

Age: _____ Weight: _____

I. Data and Calculations

A. Vital Capacity

i) **Predicted:** Use the equation below to calculate your **Predicted Vital Capacity:** _____

Equations for Predicted Vital Capacity (Kory, Hamilton, Callahan: 1960)		Where
Male	$V.C. = 0.052H - 0.022A \pm 3.60$	V.C. Vital Capacity in liters H Height in centimeters A Age in years
Female	$V.C. = 0.041H - 0.018A \pm 2.69$	

ii) **Observed:** Use the P-P result to note **Observed Vital Capacity:** _____

=

iii) Observed vs. Predicted

What is Subject's observed Vital Capacity to predicted Vital Capacity as a percentage?

Observed/Predicted VC = _____ x 100 = _____ %

Note: Vital capacities are dependent on other factors besides age and height. Therefore, 80% of predicted values are still considered "normal."

B. Volume & Capacity Measurements

Complete Table 12.2 with the requested measurement results and calculate results per the formulas provided.

Table 12.2 Measurements

Title		Measurement Result		Calculation	
Tidal Volume	TV	a = <input type="text" value="2"/> <input type="text" value="P-P"/>	Cycle 3 inhale:		$(a + b + c + d) / 4 =$
		b = <input type="text" value="2"/> <input type="text" value="P-P"/>	Cycle 3 exhale:		
		c = <input type="text" value="2"/> <input type="text" value="P-P"/>	Cycle 4 inhale:		
		d = <input type="text" value="2"/> <input type="text" value="P-P"/>	Cycle 4 exhale:		
Inspiratory Reserve Volume	IRV	<input type="text" value="2"/> <input type="text" value="Delta"/>			
Expiratory Reserve Volume	ERV	<input type="text" value="2"/> <input type="text" value="Delta"/>			
Residual Volume	RV	<input type="text" value="2"/> <input type="text" value="Min"/>			Default = 1 (Preference setting)
Inspiratory Capacity	IC	<input type="text" value="2"/> <input type="text" value="Delta"/>			TV + IRV =
Expiratory Capacity	EC	<input type="text" value="2"/> <input type="text" value="Delta"/>			TV + ERV =
Functional Residual Capacity	FRC				ERV + RV =
Total Lung Capacity	TLC	<input type="text" value="2"/> <input type="text" value="Max"/>			IRV + TV + ERV + RV =

C. Observed vs. Predicted Volumes

Using data obtained for Table 12.2, compare Subject's lung volumes with the average volumes presented in the Introduction.

Table 12.3 Average Volumes vs. Measured Volumes

Volume Title	Average Volume	Measured Volume
Tidal Volume TV	Resting subject, normal breathing: TV is approximately 500 ml. During exercise: TV can be more than 3 liters	greater than equal to less than
Inspiratory Reserve Volume IRV	Resting IRV for young adults is males = approximately 3,300 ml females = approximately 1,900 ml	greater than equal to less than
Expiratory Reserve Volume ERV	Resting ERV for young adults is males = approximately 1,000 ml females = approximately 700 ml	greater than equal to less than

II. Questions

D. Why does predicted vital capacity vary with height?

E. Explain how factors other than height might affect lung capacity.

F. How would the volume measurements change if data were collected after vigorous exercise?

G. What is the difference between volume measurements and capacities?

H. Define **Tidal Volume**.

I. Define **Inspiratory Reserve Volume**.

J. Define **Expiratory Reserve Volume**.

K. Define **Residual Volume**.

L. Define **Pulmonary Capacity**.

M. Name the **Pulmonary Capacities**.
