

## V. DATA ANALYSIS

### FAST TRACK DATA ANALYSIS

#### 1. Enter the **Review Saved Data** mode.

- Note Channel Number (CH) designations:

Channel	Displays
CH 1	Lead I
CH 2	Lead III
CH 40	Lead II (calculated)

- Note measurement box settings:


Channel	Measurement
CH 1	Delta
CH 2	Delta
CH40	Delta

#### 2. Set up the display window for optimal viewing of the first data recording.

### DETAILED EXPLANATION OF DATA ANALYSIS STEPS

If entering **Review Saved Data** mode from the Startup dialog or lessons menu, make sure to choose the correct file.

**Note:** After **Done** was pressed in the final recording section, the program used Einthoven's Law to automatically calculate Lead II from Leads I and III and added a channel for Lead II to the initial two channel recording (Fig. 6.15).



**Fig. 6.15 Example data**


The measurement boxes are above the marker region in the data window. Each measurement has three sections: channel number, measurement type, and result. The first two sections are pull-down menus that are activated when you click them.

#### Brief definition of measurements:

**Delta:** Computes the difference in amplitude between the first point and the last point of the selected area. It is particularly useful for taking ECG measurements, because the baseline does not have to be at zero to obtain accurate, quick measurements.

**Rate Mean:** If CH 40 Heart Rate data was recorded, use the Rate Mean measurement, which is designed specifically for rate data and calculates accurate statistical means using one value only for every cardiac cycle. This avoids any unintentional weighting due to time variation in heart rate, unlike the amplitude "Mean" measurement.

The selected area is the area selected by the **I-beam** tool (including endpoints).

**Note:** The append event markers  mark the beginning of each recording. Click on (activate) the event marker to display its label.



#### Useful tools for changing view:

**Display menu:** Autoscale Horizontal, Autoscale Waveforms, Zoom Back, Zoom Forward

**Scroll Bars:** Time (Horizontal); Amplitude (Vertical)

**Cursor Tools:** Zoom Tool

**Buttons:** Overlap, Split, Show Grid, Hide Grid, -, +

**Hide/Show Channel:**  + click (Windows) or  + click (Mac) the channel number box to toggle channel display.

The data window should resemble Fig. 6.16.

**Data Analysis continues...**


- Zoom in to select two consecutive clean cardiac cycles in the **Supine** recording.
- Place an **event marker** above the second R-wave to indicate which cardiac cycle will be used for measurements.

- Use the **I-Beam** cursor to select the area from the midpoint between the cycles (baseline) and the R-wave of the second cycle

 A, B

- Scroll to the **Seated** recording, select two consecutive cardiac cycles and repeat the procedure described in Steps 4 and 5.

 **Data Analysis continues...**



**Fig. 6.16 Supine data**

A clean cardiac cycle has ECG components that are easy to discern (Fig. 6.17).



**Fig. 6.17 Zoom in on Supine data**

Insert the event marker directly above the R wave of the second cardiac cycle in the display.

To place an event marker (inverted triangle,) Right-click with the cursor in the event marker region and choose the contextual menu item **Insert New Event.** You can then move the marker by holding down the **Alt** key while clicking and dragging it.

Type **Reference 1** to label the marker.

Start at the midpoint between the T-wave of cardiac cycle 1 (left) and the P-wave of cardiac cycle 2. Press and hold the mouse and sweep the cursor to the right until the end of the selected area is at peak of the desired wave— monitor the Delta measurement to determine when the actual peak is reached; small movements to the right or left may be necessary.



**Fig. 6.18 Selection from baseline to R wave peak**

**Note** R-waves may be inverted on some of the leads; include the polarity of the Delta result in the Data Report tables.

Do not use a section between the **Start of inhale** and **Start of exhale** event markers.

**Note** All remaining measurements are taken on Lead I and Lead III only so you may choose to hide Lead II (CH 40).

- Scroll to the **Start of inhale** section and select two consecutive cardiac cycles and repeat the procedure described in Steps 4 and 5.



- Scroll to the **Start of exhale** section and select two consecutive cardiac cycles and repeat the procedure described in Steps 4 and 5.



- Go back to the **Reference 1** marker created in Step 4.

- Measure the waves of the **QRS** complex and record the amplitudes for **Lead I** and **Lead III**.



- Fill in the vectorgrams.
- Answer the questions at the end of the Data Report.
- Save** or **Print** the data file.
- Quit** the program.


**END OF DATA ANALYSIS**

Type **Reference 2** to label the marker.

Type **Reference 3** to label the marker.

Use the event marker left and right arrows to move to different markers.

To measure a wave, select the area from the baseline (Isoelectric Line) to the peak of the wave.



**Fig. 6.19 Sample measurement of the Q wave**



**Fig. 6.20 Sample measurement of the R wave**



**Fig. 6.21 Sample measurement of the S wave**

- The L06 Vectorgrams are contained in the printed manual, or can be printed directly from the Help menu.
- An electronically editable **Data Report** is located in the journal (following the lesson summary,) or immediately following this Data Analysis section. Your instructor will recommend the preferred format for your lab.

**END OF LESSON 6**

Complete the Lesson 6 Data Report that follows.

## ELECTROCARDIOGRAPHY II

- *Bipolar Leads (Leads I, II, III,) Einthoven’s Law, and*
- *Mean Electrical Axis on the Frontal Plane*

### DATA REPORT

Student’s Name: \_\_\_\_\_

Lab Section: \_\_\_\_\_

Date: \_\_\_\_\_

### Subject Profile

Name: \_\_\_\_\_

Height: \_\_\_\_\_

Age: \_\_\_\_\_

Gender: Male / Female

Weight: \_\_\_\_\_

### I. Data and Plots

#### A. Einthoven’s Law—Simulated Confirmation: Lead I + Lead III = Lead II

Table 6.1 Supine

Lead	Same Single Cardiac Cycle	mV*
Lead I	1 Delta	
Lead III	2 Delta	
Lead II	40 Delta	

\*Include the polarity (+ or -) of the Delta result since R-waves may be inverted on some of the leads.

#### B. Mean Electrical Axis of the Ventricles (QRS Axis) and Mean Ventricular Potential—Graphical Estimate

Use Table 6.2 to record measurements from the Data Analysis section:

Table 6.2


CONDITION	QRS	
	Lead I 1 Delta	Lead III 2 Delta
Supine		
Seated		
Start of inhale		
Start of exhale		

One way to approximate the mean electrical axis in the frontal plane is to plot the magnitude of the R wave from Lead I and Lead III, as shown in the Introduction (Fig. 6.4).

1. Draw a perpendicular line from the ends of the vectors (right angles to the axis of the Lead) using a protractor or right angle guide.
2. Determine the point of intersection of these two perpendicular lines.
3. Draw a new vector from point 0.0 to the point of intersection.

The direction of this resulting vector approximates the mean electrical axis (QRS Axis) of the ventricles. The length of this vector approximates the mean ventricular potential.

Create two plots on each of the following graphs, using data from Table 6.2. Use a different color pencil or pen for each plot.



From the above graph, find the following values:

Condition	Mean Ventricular Potential	Mean Ventricular (QRS) Axis
Supine	_____	_____
Seated	_____	_____

Explain the difference (if any) in Mean Ventricular Potential and Axis under the two conditions:

---




---



---



---



From the above graph, find the following values:

Condition	Mean Ventricular Potential	Mean Ventricular (QRS) Axis
Start of inhale	_____	_____
Start of exhale	_____	_____

Explain the difference (if any) in Mean Ventricular Potential and Axis under the two conditions:

---



---



---



---


**C. Mean Electrical Axis of the Ventricles (QRS Axis) and Mean Ventricular Potential—More Accurate Approximation**

Use Table 6.3 to add the Q, R, and S potentials to obtain net potentials for Recording 16 Supine.

Table 6.3

POTENTIAL	QRS	
	Lead I <input type="text" value="1"/> <input type="text" value="Delta"/>	Lead III <input type="text" value="2"/> <input type="text" value="Delta"/>
Q		
R		
S		
QRS Net		

Graph 3: *Supine*



From the above graph, find the following values:

Condition	Mean Ventricular Potential	Mean Ventricular (QRS) Axis
Supine	_____	_____

Explain the difference in Mean Ventricular Potential and Axis for the Supine data in this plot (Graph 3) and the first plot (Graph 1).

---



---



---

**II. Questions**

D. Define ECG.

---

---

---

---

E. Define Einthoven’s Law.

---

---

---

---

F. Define Einthoven’s Triangle and give an example of its application.

---

---

---

---

G. What normal factors effect a change the orientation of the Mean Ventricular (QRS) Axis?

---

---

---

---

H. Define Left Axis Deviation (LAD) and its causes.

---

---

---

---

I. Define Right Axis Deviation (RAD) and its causes.

---

---

---

---

J. What factors affect the amplitude of the R wave recorded on the different leads?

---

---

---

---



**III. OPTIONAL Active Learning Portion**

A. *Hypothesis*

---

---

---

---

B. *Materials*

---

---

---

---

---

---

C. *Method*

---

---

---

---

---

---

D. *Set Up*

---

---

---

---

---

---

E. *Experimental Results*

---

---

---

---

---

---