

AFT13 Disposable Pulmonary Function Filter and Mouthpiece *Available in packs of 10 or 250*

NOTE: The AFT13 was discontinued in June of 2016. For SS11L/LA, see AFT1filter + AFT2 mouthpiece. For SS11LB, see AFT36 for current product offering.

Eliminate cross-contamination concerns with this bacteriological filter with disposable plastic-coated paper mouthpiece to protect subjects and equipment. These exceed all recommended performance standards with 99.9% bacterial filtration efficiency and 99.9% viral filtration efficiency. They feature low resistance and minimal dead space (45 ml when measured without tube fittings). These surpass published ATS recommendations for flow resistance in pulmonary function instrumentation, which suggest resistance should be below 1.5 cm H₂O/L/sec at flow rates less than 12 L/sec. Port: 30 mm OD.

Requires AFT36 + AFT11D or AFT4 + AFT8 + AFT11D for proper operation.