

AFT26 Calibration Syringe (2.0 liter)

NOTE: The AFT26 2.0 liter calibration syringe was discontinued in September of 2017. For current offering, see the AFT27 3.0 liter syringe.



The AFT26 is a 2.0 Liter Calibration Syringe for the SS11LB, SS11LA or TSD117A Airflow Transducer. The AFT26 Calibration Syringe is certified to have a 2-liter volume that meets or exceeds an accuracy $\pm 1\%$ of the total displacement volume. The increased size and accuracy of this 2.0 liter calibration syringe provide a wider calibration range than the AFT6A for advanced studies. A coupler is included and can be

reordered as AFT11D (SS11LB) or AFT11I (SS11LA) if it is inadvertently discarded when an airflow accessory is removed.

AFT11I COUPLER

NOTE: The AFT11I coupler for connecting the AFT26 2.0 Liter Calibration Syringe to the SS11LA airflow transducer was discontinued in September of 2019. To connect an existing AFT26 syringe to the SS11LA, use the following connections:

AFT26 + AFT11B + AFT11D + SS11LA