

CBLEPM**3.5 MM TO 2 X TINNED WIRE**

The CBLEPM connects the following:

- STMISOLA/L to a Measurement Computing 4-channel D/A unit. 3.5 mm phone jack connects to STMISOLA/L phone plug trigger and two tinned wires connect to the D/A unit screw terminals.
- FOTS100A module to BIOPAC Research System (additional cables also required, see module page for details)
- MPMS200A module to BIOPAC Research System (additional cables also required, see module page for details)

**SPECIFICATIONS**

3.5 mm phone plug (female) to 2 x tinned wires

See also: [STMISOLA Stimulator](#) and [STMEPM E-Prime System](#)