




Inspiring People and Enabling Discovery About Life

42 Aero Camino  
Goleta, CA 93117 USA  
T: 805.685.0066  
F: 805.685.0067  
info@biopac.com  
www.biopac.com

## WORKAROUND FOR APPLE SECURITY ERROR


### The Error:

Apple cannot check the installer file for malicious software.




### The Workaround:

Click on the Apple icon in the top left of the desktop, go into *System Preferences*, and choose *Security & Privacy* (make sure to unlock to initiate changes, at the bottom left of the menu there should be a padlock icon).



## WORKAROUND FOR APPLE SECURITY ERROR

Under the *General* tab, at the very bottom, you will see the dialog “Allow apps downloaded from.”



Choose the option “App Store and identified developers.”

A small line of text should display below saying it still does not recognize the software as an “identified developer.” Click “Open Anyway”.

- If this additional text is not there, try opening the .pkg file again.

At this point the computer will run through a verification process. A similar message with the malicious software warning may appear again, but this time it will give the option to open the software (rather than just “OK”). Click “Open.”

