

### ETHUSB USB 2.0 ETHERNET ADAPTER



*ETHUSB is included in MP160/150 Systems, upgrades to MP160/150, and VR Systems.*

Use to connect to a 10/100 Mbps network through a USB port—no need to open up your computer case to add an internal Ethernet card. The adapter's compliance with USB 2.0 (480 Mbps) ensures true 10/100 Mbps network speed without any compromise. Adapter is compact and USB bus-powered; no external power adapter required.

#### Key Features

- Instantly connect to a 10/100 Mbps network through a USB port—no need to open up your computer case to add an internal Ethernet card
- Compliant with USB 2.0 and USB 1.1 specifications
- Compliant with IEEE 802.3 (10Base-T) and 802.3u (100Base-TX) standards
- Powered by USB port—no external power adapter required
- Supports both full-duplex and half-duplex operations
- Supports suspend mode and remote wakeup via link-up and magic packet
- Equipped with diagnostic LEDs

#### System Requirements

- IBM compatible Pentium-233 MHz or faster PC or Mac
- 64 MB RAM or more
- One available USB port
- Windows 10, 8, 7, Mac OS X 10.x

#### Package Includes

- USB 2.0 to 10/100 Ethernet Adapter
- CD (Driver & User Manual)
- Quick Install Guide

#### Specifications

Cable Type: USB  
Dimensions: 6.60 cm x 2.27 cm x 1.52 cm [ 2.60" x 0.90" x 0.60" ]  
Weight: 0.03 kg [ 0.06 lbs. ]

Connector A: USB A (male)



Connector B: RJ45 (female)



Certifications: 802.3; 802.3u; USB 1.1; USB 2.0