

ETHERNET ACCESSORIES**CBLETH1/2 – Ethernet Cables**

CBLETH1 is a 2-meter Ethernet patch cable.

- Use one CBLETH1 to connect the MP150 to an Ethernet Switch and one CBLETH1 to connect the Switch to a local area network (LAN).



CBLETH2 is a 2-meter Ethernet crossover cable.

- MP150 Systems include one CBLETH2.
- Use the CBLETH2 to connect the MP150 to an Ethernet interface (such as the ETHUSB interface shipped with your MP System).

