

EXTENSION CABLES**CBLEXT** Serial Extension Cable

Use this 3.6-meter extension cable to increase the distance between the MP100 or MP30 and your computer. The CBLEXT is used to extend the length of CBLSERA (use only one CBLEXT per MP100 or MP30 System).

**CBLHLT1** Phono Plug Extension Cable

Use this 7.6-meter extension cable with phono connectors for

- AMI100D/HLT100C high level transducer module.
- NIBP100D Noninvasive Blood Pressure System (cable included) between the TCI105 and NIBP100D connector.

**CBLOXY-EXT** Pulse Oximeter Extension Cable

Use this 2.4-meter cable to extend the connection from an OXY100E or OXY200 SpO₂ amplifier module to a TSD124 series human or TSD270 series veterinary SpO₂ transducer. (Replaces the OXY100E-200-EXT cable)

