



FACEREADER-ACK – FACEREADER WITH ACQKNOWLEDGE

Real-time expression data synchronized with Physiology data and/or eye tracking



FaceReader™ is the premier professional software for automatic analysis of facial expressions. Emotion reading software locks onto a subject's face and analyzes facial movements to classify subject's response. Easily integrate with physiological parameters from an MP160 Research System or eye tracking data. Start *AcqKnowledge* to begin recording, and FaceReader data is automatically synchronized and recorded in the same graph file. Monitor data in the FaceReader display for real-time feedback. Classifications include:

- Happy
- Sad
- Scared
- Disgusted
- Surprised
- Angry
- Contempt
- Neutral



Seamless Integration – easily integrate and synchronize physiological data or eye tracking through *AcqKnowledge*.

Complete Facial Expression Analysis – accurate modeling of the face by describing over 500 key points delivers objectivity in observations.

Useful for a Variety of Applications – psychology, educational research, consumer behavior research, usability studies, Neuromarketing research, and more!

System Requirements

Microsoft® Windows 10, 8, 7, 64-bit editions, professional workstation with high-end video supporting Direct3D acceleration.

[Click to request more info about FaceReader.](#)