

RX137 SERIES REPLACEMENT AIRFLOW HEADS**For TSD137 Series Pneumotachs**

The RX137 series are low airflow heads for the TSD137 series pneumotach transducers. The RX137 heads can be mixed and matched with any of the TSD137 series pneumotachs. Switching one head for another when using a single TSD137 pneumotach can accommodate a wide range in flows. RX137 heads connect to the TSD160A differential pressure transducer via standard 3mm or 4mm ID tubing. Multiple RX137 heads help eliminate equipment downtime during cleaning procedures.

RX137 SERIES SPECIFICATIONS

Head	Range (ml/sec)	Dead Space (cc)	Length (mm)	Flow Ports		Weight (grams)
				ID (mm)	OD	
RX137A	±12	0.1	75	1.35	7	100
RX137B	±20	0.8	75	6.00	7	90
RX137C	±60	0.9	75	6.00	7	90
RX137D	±150	2.0	75	9.00	10	100
RX137E	±350	4.0	60	10.00	11	60