

ELECTRODE LEAD SETS SS2L and SS2LA

- **These items are discontinued—see SS2LB for current product**

General-purpose electrode cables for almost all applications requiring the use of electrodes. One end of the SS2L/SS2LA cable has a Smart Sensor connector on it that connects to an MP3X/45 data acquisition unit and the other end splits into three smaller cables; each end of the smaller cables is fitted with a connector that clamps onto electrodes.

Used to connect disposable electrodes placed on the surface of the skin to MP3X/4X data acquisition unit to measure—depending on electrode placement—muscle contraction, heartbeats, brainwaves, or other biopotentials.

SS2L and SS2LA Specifications

Cable Length: 2 meters
Connector Type: pinch connector
Interface connector: 9 Pin DIN



use

an