

SS70L Isolated BNC Input Adapter for MP36/MP35

NOTE: The SS70L Isolated BNC Input Adapter was discontinued in June of 2017. See [SS70LA](#) for current product offering.



This BNC adapter is required when connecting un-isolated third party devices (i.e. amplifiers, chart recorders or signal generators), while electrodes, attached to human Subjects are connected to other input channels.

Connector Type: BNC

Signal range: ± 20 V (MP36/MP36R/MP45)
 ± 10 V (MP35)

WARNING! Since all MP inputs share a common isolated ground, connecting an un-isolated device to any channel voids the isolation for all channels and exposes the Subject to possible shock hazards.