


**M-POD EVENT PRESENTATION & MARKING SYSTEMS**

STM-M-POD-SYS-IO (for MP36/36R Systems)

STM-M-POD-SYS-STP (for MP160 Systems)



These Event Presentation & Marking Systems include m-pod, SuperLab Pro Stimulus Presentation Software, Response Pad (RB-740), and BIOPAC interface cable (option to interface MP160 System also includes Isolated Digital Interface STP100C ).

SuperLab offers a host of powerful features, including

- Playing movies
- Stimulus lists
- Support for JPEG, GIF, PNG, and TIFF files
- Built-in support for RSVP and self-paced reading
- Improved support for fMRI and EEG/ERP
- Trial variables
- Conditional branching (if/then/else)
- Multiple input devices in the same experiment
- Unicode application that handles Japanese, Chinese, and other international fonts just as easily as it handles English fonts.

**See also:** Product description and specs for m-pod, StimTracker, SuperLab, and STP100C.