


BIOHARNESS TEAM SYSTEMS

TEAMSYSTEM-4A – Team System Physiology Monitoring for 4 Subjects



TEAMSYSTEM-4A

TEAMSystem-4A is an RF telemetry physiology monitoring system with the same functionality as TEAMSystem-1A but expanded for 4 people. Includes x 4 BioHarness units, 4-unit charger and any size combination of 4 straps.

TEAMSYSTEM-4A Components

- BioHarness Unit x 4 for TEAMSystem-4A
- BioHarness Charger (holds 4 units in TEAMSystem-4A)
- Smart Fabric™ strap - specify size when ordering: Small fits 71-97 cm (28-38 inches) or Large fits 91-122 cm (36-48 inches) (TEAMSystem-4A – specify any combination of four straps)
- USB Power
- 5 m USB ext cable
- 1 m USB A to miniB cable
- Support Docs (Quick Guide, Install Guide, User Guide)
- Team System CD
- System Case
- Antenna with stand and case

The TEAMSystem-4A was discontinued in April of 2016. For current solutions, see [TEAMSystem-10](#), [TEAMSystem-30](#), or [TEAMSystem-50](#).