


TSD124 SERIES SPO2 TRANSDUCERS FOR OXY100E



TSD124 Series



TSD124A Finger



TSD124C Flex Wrap

The TSD124 series human oximetry transducers are reliable and simple to use on a wide range of subjects for both short-term and continuous noninvasive monitoring. The transducers incorporate Nonin's PureLight® sensors and are backed by a six-month warranty. Use with the OXY100E oximetry amplifier.

AVAILABLE TYPES:

TSD124A Finger Clip Transducer - The TSD124A has been discontinued and replaced by the TSD124D.

Subject Range: > 30 kg (66 lbs)

Preferred application: Index, middle or ring fingers

TSD124B Ear Clip SpO₂ Transducer

Subject Range: > 40 kg (88 lbs)

Length: 1 m

TSD124C Flex Wrap SpO₂ Transducer (Ships with 25 adhesive wrap guides)

Also available: OXY100E-200 EXT Pulse Oximeter extension cable – 3m