

RX TISSUE BATH ACCESSORIES / REORDER PARTS


Tissue Holders



Tissue Clips



Warming Coil



Oxygen Filter



Tissue Bath




Reservoir



Mount Accessories



Field Stimulation Electrode



RXHOLDER-S

Tissue Holder (stainless steel)

RXHOLDER-G

Tissue Holder (glass)

RXHOLDER-T

Triangle Tissue Holder (stainless)

RXCLIP

Tissue Clip (stainless steel)

RXCLIP-TRI

Triangle Tissue Clip for Rings (stainless steel)

RXCOIL

Warming Coil

RXO2FILTER

Oxygen Filter (glass)

RXBATH

Tissue Bath

RXRESERV

Reservoir 1000ml

RXMOUNT

Mount Accessories Kit

STIMHOLDER

Field Stimulation Electrode

TISSUE BATH ACCESSORIES SPECIFICATIONS

1 x Tissue Holder—stainless steel; 15 mm high x 9 mm wide; reorder as RXHOLDER-TR

1 x Tissue Holder—glass; reorder as RXHOLDER-G

1 x Tissue Holder —stainless steel; reorder as RXHOLDER-S

2 x Tissue Clip—stainless steel; 15 mm high x 5 mm wide: reorder as RXCLIP

2 x Triangle Tissue Clip—stainless steel; 15 mm high x 12 mm wide; reorder as RXCLIP-TRI

1 x Integrated heater—1,600 ml volume, programmable temp. 20° - 44° C

1 x Circulator pump—15 W; 500 ml/min

1 x Warming Coil; reorder as RXCOIL

1 x Oxygen Filter; reorder as RXO2FILTER

1 x Bath —reorder as RXBATH

1 x Reservoir—1000 ml; reorder as RXRESERV

Mount Accessories Kit; reorder as RXMOUNT

Field Stimulation Electrode; reorder as STIMHOLDER

1 x Micrometer-transducer assembly

1 x 3 way Rotary Valve

1 x Power Supply – 110V/60 Hz or 220V/50 Hz